

IBUMAMZ-

The unique mounting system of the sidewalls, positioned inside the perimeter of the filter, prevents contaminants from entering into the work area through the same filter frame thanks to the Venturi-effect.

The air passes through a polyurethane fiber G-4 prefilter and it is pushed by the fan to the retention plenum, then goes through the absolute H-14 HEPA filter and reaches the working chamber. The airflow leaving the H-14 filter placed in front of the operator, at a constant speed, creates an ISO 5 sterile area removing any possible contaminants from the product you are working with.

USES

- \checkmark Parenteral bags reconstruction
- \checkmark Assisted reproductive technologies
- √ Micropropagation
- √ Plant Cell Cultures
- √ Food and water quality control
- \checkmark Dust Protection in optics and electronics industry

NEW TECHNICAL FEATURES



More info on the new LCD screen

- √ Large format 127x34mm
- \checkmark Type of filter installed, working hours, expiry date and next revision
- √ Countdown Timer
- \checkmark Clock and calendar

New features and components

- \checkmark Initial cycle flow adequacy and final cycle purge
- √ Fault LED indicator
- √ Air Flow Control by Microprocessor
- $\checkmark\,$ HEPA filters with electronic chip
- $\checkmark\,$ LED Lighting

New alarms and warnings

- √ Next Validation
- \checkmark Few hours of filter life
- √ Countdown timer
- $\sqrt{}$ Filter has expired (hours of use)
- $\sqrt{}$ Filters has expired (expiry date)
- √ Cabinet without filter

It is not a typographical error, *7 year warranty*

Because we are convinced of the quality of our products.



More information www.cruma.es



Illumination level Noise Level Packaging: 10 with internation Well done, well shipped. Our responsible packaging

Wood box 100% recyclable with international phytosanitary <u>c</u>ertificate.



www.cruma.es export@cruma.es



CRUMAHZ-2

✓ Manufacture: outer surfaces made of cold rolled steel gray / ocher coating, stoveenamelled, resistant to the chemical compounds used for sterilization

2 20000


- $\checkmark\,$ Worktop: the surface of the working chamber is in stainless steel AISI 304 2B
- √ Acrylic transparent side walls
- ✓ Absolute HEPA H-14 filter with flow equalizer, efficiency greater than 99.995% on 0.3 micron particles. Filter is in accordance with EN-779:1996 and EN-1822:1998
- $\sqrt{\rm \ Fan:}$ centrifugal double suction, suitable to provide laminar flow in the workspace
- ✓ Control panel with soft-touch keyboard and switch for LED and optional UV lamps; the air speed is kept at a constant rate thanks to the HEPA filter which is electronically self- regulated, electronic control air speed and control key to slow it down
- √ Filter life timer
- √ Lighting: LED lamp placed in the central part of the cabinet, for a better light distribution over the working surface

TECHNICAL INFORMATION Number of filtering columns 1 2 Number of IP44 fans 0,50 m/s Average face velocity Total electrical power consumption 406 W Voltage/Frequency 110-220 V / 50-60 Hz Illumination level LED 18 W / 900 Lux 45 dB Noise Level Volume 0,65 m³ Packaging: 100% recyclable Wood box with international phytosanitary certificate Weight 120 Kg




Contact your distributor or call us if you have any questions or need technical support, spare parts, maintenance service... +34 93 370 61 62







SIZES (mm)



SIZES (mm)						
Exterior			Interior			
Width 1295	Depth 770	Height 977	Width 1169	Depth 600	Height 633	

STRUCTURE	
Metallic parts	1.2 mm galvanized steel, coated with antiacid polymer powder resin thermo-hardened at 200 ° C
Side Panels	Transparent polymethylmethacrylate 8 mm thick (light transmission of 93%)





HORIZONTAL LAMINAR FLOW CABINET TO PROTECT THE SAMPLE





STANDARD EQUIPMEN	Т
Electronic circuit with large format LCD display	Security levels: level 1 for users and level 2 for maintenance users
Filters controlled by microprocessor	The filters are provided with a microchip with mini USB connection identifying the type of filter installed, the expiry date and the serial number
Illumination	96 LED high light intensity and low consumption 18 Watts / 900 Lux
Digital timer with audible warning	Hourmeter for the control and warning of the work in the cabinet
Electrical plug	Placed in the internal roof of the cabinet and conmuted with the electronic circuit
Clock and Calendar	Display of date and time
Steel work surface	Worktop of stainless steel surface AISI 304 2B
G4 prefilter	Pre-filtering class G4 biofibre synthetic blanket
Warranty	7 years



OPTIONAL EQUIPMENT

Tubular stand

Portable UV light (2 x 15 Watt UV lamps)

COMPLIANCE TO STANDARDS

Filters	EN-779: 1996 (HEPA & ULPA Filters) EN-1822:1998 (HEPA & ULPA Filters)
Quality	UNE EN ISO 9001:2015

We recognise our responsability and dependence towards a healthy environment and, therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

our 🕫 engagement





