



PRODUCT INFORMATION

Product VECTASHIELD® Mounting Medium

Catalog No. H-1000

Amount 10 ml

Storage Conditions Refrigerate in the dark

VECTASHIELD® Mounting Medium is a unique, stable formula for preserving fluorescence. VECTASHIELD® Mounting Medium prevents rapid photobleaching of fluorescent proteins and fluorescent dyes (see website for compatibility).

VECTASHIELD® Mounting Medium is supplied in a bottle fitted with a screw cap. A cap fitted with a drop dispenser pipet is also included. The screw cap can be replaced with the drop dispenser cap for routine use and during storage. Once the drop dispenser cap has been secured on the bottle, the bottle should be stored in an upright position.

To mount tissues or cells on a slide, dispense VECTASHIELD® Mounting Medium onto the specimen. The pipet supplied with the product is designed to deliver small drop volumes of approximately 25 µl. Coverslip and allow VECTASHIELD® Mounting Medium to disperse over the entire section.

VECTASHIELD® Mounting Medium does not solidify, but remains a liquid on the slide and can be stored without sealing. If desired, coverslips can be sealed around the perimeter with nail polish or a plastic sealant. Mounted slides should be stored at 4 °C, protected from light.

VECTASHIELD® Mounting Media Products:

Product	Counterstain	Cat. No.	Unit Size	Hardening	Refractive Index*
VECTASHIELD® Mounting Medium	none	H-1000	10 ml	no	1.45
	DAPI	H-1200	10 ml	no	1.45
	PI	H-1300	10 ml	no	1.45
VECTASHIELD® Hard•Set™ Mounting Medium	none	H-1400	10 ml	yes	1.46
	DAPI	H-1500	10 ml	yes	1.46
	TRITC-Phalloidin	H-1600	10 ml	yes	1.46

*Measured at room temperature. VECTASHIELD® Hard•Set™ Mounting Medium measured after hardening.

VECTASHIELD is a registered trademark of Vector Laboratories, Inc.

01/15