



**LECTINS
AND GLYCOBIOLOGY**

STATE OF THE ART LABELING & DETECTION SYSTEMS



Vector Laboratories Inc., USA
30 Ingold Road
Burlingame, CA 94010
(800) 227-6666 (orders)
(650) 697-3600 (technical service)
(650) 697-0339 (fax)
vector@vectorlabs.com (email)
www.vectorlabs.com

Vector Laboratories Ltd., UK
3 Accent Park, Bakewell Road
Orton Southgate, Peterborough PE2 6XS
(01733) 237999 (telephone)
(01733) 237119 (fax)
vector@vectorlabs.co.uk (email)
www.vectorlabs.com


Vector Laboratories, Inc., Canada
3390 South Service Road
Burlington, Ontario L7N 3J5
(888) 629-2121 (orders only)
(905) 681-0900 (tel)
(905) 681-0166 (fax)
vectorca@netinc.ca (email)
www.vectorlabs.com

Introduction

Since the 1880's, it has been known that extracts from certain plants could agglutinate red blood cells. In the 1940's, agglutinins were discovered which could "select" types of cells based on their blood group activities. Although "lectin" was originally coined to define agglutinins that could discriminate among types of red blood cells, today the term is used more generally and includes sugar-binding proteins from many sources regardless of their ability to agglutinate cells. Lectins have been found in plants, viruses, microorganisms, and animals but despite their ubiquity, in many cases their biological function is unclear.

Most lectins are multimeric, consisting of non-covalently associated subunits. It is this multimeric structure that gives lectins their ability to agglutinate cells or form precipitates with glycoconjugates in a manner similar to antigen-antibody interactions. This unique group of proteins has provided researchers with powerful tools to explore a myriad of biological structures and processes. Because of the specificity that each lectin has toward a particular carbohydrate structure, even oligosaccharides with identical sugar compositions can be distinguished or separated. The affinity between a lectin and its receptor may vary a great deal due to small changes in the carbohydrate structure of the receptor. These properties enable the researcher to discriminate between structures, isolate a specific glycoconjugate, cell, or virus from a mixture, or study one process among several. Another property of some lectins is an ability to induce mitosis in cells that are normally not dividing. This property has been exploited extensively in an attempt to understand the process of lymphocyte blastogenesis and the biochemical and structural alterations associated with mitogenesis. These important applications include lymphokine production and viral coat protein (eg gp120) isolation. In fact, thousands of articles on lectins have been published examining hundreds of different aspects and uses of lectins.

Lectins have been purified by "conventional" procedures including salt-induced crystallization, ethanol precipitation, ion exchange chromatography and gel filtration, or by affinity chromatography. The former methods rely on the physicochemical properties of the proteins for separation while affinity chromatography depends on the specific interaction between the lectin and a carbohydrate structure attached to an inert matrix. We employ both "conventional" procedures and affinity chromatography for each of our lectins. Purification is monitored and final product is assessed by immunoprecipitation with antisera, agglutination titre, polyacrylamide gel electrophoresis, and binding activity to specific affinity columns, providing the assurance that our customers have the best lectins available.



Mouse Tongue: endothelial cells stained with Dylight® 594-labeled Griffonia simplicifolia Lectin, Isolectin B₄ (red fluorescence). Mounted with VECTASHIELD® Hard•Set™ with DAPI (blue fluorescence).

Table of Lectin Properties

Lectin	Common Abbreviation	Source	Mol. Wt. (kDa)	Number of Subunits	pI	ϵ 0.1% 280 nm	Glycoprotein	Metal Ions Present
<i>Agaricus bisporus</i>	ABL	<i>Agaricus bisporus</i> white button mushrooms	68	4	5-6	1.6	–	No
<i>Aleuria aurantia</i>	AAL	<i>Aleuria aurantia</i> mushrooms	72	2	9	2.97	No	--
<i>Amaranthus caudatus</i>	ACL, ACA	<i>Amaranthus caudatus</i> seeds	66-70	2	6.7 - 7.7	1.60	No	No
<i>Bauhinia purpurea</i>	BPL, BPA	<i>Bauhinia purpurea alba</i> (Camel's Foot Tree) seeds	195	4	4.6 - 6	1.75	Yes	No
Concanavalin A	Con A	<i>Canavalia ensiformis</i> (Jack Bean) seeds	104	4	6.3 - 7	1.20	No	Ca ⁺⁺ , Mn ⁺⁺
Succinylated Concanavalin A	Succinylated Con A	<i>Canavalia ensiformis</i> (Jack Bean) seeds	56	2	< 4.4	1.2	No	Ca ⁺⁺ , Mn ⁺⁺
<i>Datura stramonium</i>	DSL	<i>Datura stramonium</i> (Thorn Apple, Jimson Weed) seeds	86	1	> 9	0.80	Yes	No
<i>Dolichos biflorus</i>	DBA	<i>Dolichos biflorus</i> (Horse Gram) seeds	111	4	4.6 - 5	1.22	Yes	Ca ⁺⁺ , Mn ⁺⁺ , Mg ⁺⁺ , Zn ⁺⁺
<i>Erythrina cristagalli</i>	ECL, ECA	<i>Erythrina cristagalli</i> (Coral Tree) seeds	54	2	6.3 - 6.5	1.30	Yes	Ca ⁺⁺ , Mn ⁺⁺ , Zn ⁺⁺
<i>Euonymus europaeus</i>	EEL	<i>Eunonymus europaeus</i> (Spindle Tree) seeds	140	4	4.4	2.40	Yes	Ca ⁺⁺ , Zn ⁺⁺
<i>Galanthus nivalis</i>	GNL	<i>Galanthus nivalis</i> (Snowdrop) bulbs	50	4	3.5 - 4	1.90	No	No
<i>Griffonia (Bandeiraea) simplicifolia</i> I	GSL I, BSL I	<i>Griffonia (Bandeiraea) simplicifolia</i> seeds	114	4	5 - 6.5	1.40	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Griffonia (Bandeiraea) simplicifolia</i> I Isolectin B ₄	GSL I - B ₄	<i>Griffonia (Bandeiraea) simplicifolia</i> seeds	114	4	6 - 6.2	1.4	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Griffonia (Bandeiraea) simplicifolia</i> II	GSL II, BSL II	<i>Griffonia (Bandeiraea) simplicifolia</i> seeds	113	2	5 - 6	1.25	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Hippeastrum hybrid</i>	HHL, AL	<i>Hippeastrum hybrid</i> (Amaryllis) bulbs	50	4	4.7 - 5.1	1.85	No	No
Jacalin	Jacalin	<i>Artocarpus integrifolia</i> (Jackfruit) seeds	66	4	7.8	1.50	Yes	No
<i>Lens culinaris</i>	LCA, LcH	<i>Lens culinaris</i> (lentil) seeds	50	4	7.6 - 8.4	1.25	No	Ca ⁺⁺ , Mn ⁺⁺
<i>Lotus tetragonolobus</i>	LTL	<i>Lotus tetragonolobus</i> , <i>Tetragonolobus purpurea</i> (Winged Pea, Asparagus Pea) seeds	107	4	7.3, 7.6, 7.9 & 8.2	1.51	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Lycopersicon esculentum</i>	LEL, TL	<i>Lycopersicon esculentum</i> (tomato) fruit	71	1	>9	0.76	Yes	--
<i>Maackia amurensis</i> I	MAL I, MAL	<i>Maackia amurensis</i> seeds	130	2	4.7	1.40	Yes	No
<i>Maackia amurensis</i> II	MAL II, MAH	<i>Maackia amurensis</i> seeds	130	2	4.7	1.33	Yes	No
<i>Maclura pomifera</i>	MPL	<i>Maclura pomifera</i> (Osage Orange) seeds	44	4	4.8 - 5.3	1.40	No	No
<i>Musa paradisiaca</i>	BanLec	<i>Musa paradisiaca</i> (banana fruit)	30	2	7.3	0.8	–	–

Lectin	Mitogenic Activity	Blood Group Specificity	Preferred Sugar Specificity	Inhibitor or Eluting Sugar	Special Application
<i>Agaricus bisporus</i>	No	Non-specific	Gal(β-1,3) GalNAc	Fetuin	Binds sialylated and non-sialylated T antigen
<i>Aleuria aurantia</i>	No	Non-specific	Fucα6GlcNAc	L-Fuc	Tumor cell marker detection
<i>Amaranthus caudatus</i>	Yes	O>B>A (-SA) T antigen	Galβ3GalNAc	desialylated fetuin	T-antigen probe
<i>Bauhinia purpurea</i>	Yes	A,B,O (-SA)	Galβ3GalNAc	Lactose	Detection of Reed Sternberg cells of Hodgkin's disease; binds N-Antigen of human erythrocytes
Concanavalin A	Yes	Non-specific	αMan, αGlc	MeαMan+ MeαGlc	Polylysine and histone conjugated as vectors to transfer genes to airway epithelial cells; insulin-like activity
Succinylated Concanavalin A	Yes	None	αMan, αGlc	MeαMan+ MeαGlc	Growth inhibitor of 3T3 mouse fibroblasts
<i>Datura stramonium</i>	Yes	A, B, O	(GlcNAc) _{2,4}	Chitin hydrolysate	Marker for mouse peritoneal cells
<i>Dolichos biflorus</i>	No	A ₁ >>A ₂	αGalNAc	GalNAc	Distinguishes A ₁ from A ₂ human red blood cells; binds to rat macrophages in lung tissues
<i>Erythrina cristagalli</i>	Yes	A (-SA)	Galβ4GlcNAc	Lactose	Isolation of human natural killer cells
<i>Euonymus europaeus</i>	Yes	O (-SA), B	Galα3Gal	Lactose	Endothelial cell binding
<i>Galanthus nivalis</i>	No	Rabbit	αMan	MeαMan	Detection of HIV and SIV glycoprotein; isolation of mouse IgM; virus isolation
<i>Griffonia (Bandeiraea) simplicifolia I</i>	No	B>>A1	αGal, αGalNAc	Gal/GalNAc	Endothelial cell marker (for mouse)
<i>Griffonia (Bandeiraea) simplicifolia I Isolectin B₄</i>	No	B	αGal	Gal or Raffinose	Non-primate endothelial cell marker; neuronal marker
<i>Griffonia (Bandeiraea) simplicifolia II</i>	No	A (-SA)>>B (-SA)	α or βGlcNAc	Chitin hydrolysate or GlcNAc	Selective staining of Golgi apparatus; marker for uterine blood vessels and certain carcinomas
<i>Hippeastrum hybrid</i>	No	Rabbit	αMan	MeαMan	
Jacalin	Yes	O (+SA), T antigen	Galβ3GalNAc	Gal or Melibiose	Purification of human IgA
<i>Lens culinaris</i>	Yes	Non-specific	αMan, αGlc	MeαMan+ MeαGlc	
<i>Lotus tetragonolobus</i>	No	O<A2	αFuc	L-Fuc	Distinguishes between pathogenic and non-pathogenic trypanosomes
<i>Lycopersicon esculentum</i>	No	Non-specific	(GlcNAc) _{2,4}	Chitin hydrolysate	Perfusion studies in mouse (binding to vascular endothelium)
<i>Maackia amurensis I</i>	Yes	Non-specific	Galβ4GlcNAc	Lactose	
<i>Maackia amurensis II</i>	Yes	Non-specific	Neu5Acα3Galβ3GalNAc	Human Glycophorin	
<i>Maclura pomifera</i>	Yes	A, B, O (-SA)	Galβ3GalNAc	Gal	Binding to rat lymphoid cells
<i>Musa paradisiaca</i>	Yes	Rabbit	αMan, αGlc	MeαMan	Binds gp 120 glycoprotein

Sugar Abbreviations:

Fuc	L-Fucose	Man	Mannose
Gal	D-Galactose	MeαGlc	α-Methylglucoside
GalNAc	N-Acetylgalactosamine	MeαMan	α-Methylmannoside
Glc	D-Glucose	Neu5Ac	N-Acetylneuraminic acid (sialic acid)
GlcNAc	N-Acetylglucosamine	SA	Sialic Acid

Table of Lectin Properties (continued)

Lectin	Common Abbreviation	Source	Mol. Wt. (kDa)	Number of Subunits	pI	$\epsilon_{0.1\%}$ 280 nm	Glycoprotein	Metal Ions Present
<i>Narcissus pseudonarcissus</i>	NPL, NPA, DL	<i>Narcissus pseudonarcissus</i> (Daffodil) bulbs	59	4	4.2 - 4.6	1.87	No	No
Peanut	PNA	<i>Arachis hypogaea</i> peanuts	110	4	5.5 - 6.5	0.89	No	Ca ⁺⁺ , Mg ⁺⁺
<i>Phaseolus vulgaris</i> Erythroagglutinin (PHA-E)	PHA-E	<i>Phaseolus vulgaris</i> (Red Kidney Bean) seeds	126	4	6 - 8	1.16	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Phaseolus vulgaris</i> Agglutinin (PHA-E+L)	PHA-E+L	<i>Phaseolus vulgaris</i> (Red Kidney Bean) seeds	126	4	5.2 - 6.2	1.16	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Phaseolus vulgaris</i> Leucoagglutinin (PHA-L)	PHA-L	<i>Phaseolus vulgaris</i> (Red Kidney Bean) seeds	126	4	4.2 - 4.8	1.16	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Pisum sativum</i>	PSA	<i>Pisum sativum</i> (Pea) seeds	53	4	6.0 - 6.7	1.20	Trace	Ca ⁺⁺ , Mn ⁺⁺
<i>Psophocarpus tetragonolobus</i> I	PTL I, WBA I	<i>Psophocarpus tetragonolobus</i> (Winged Bean) seeds	57	2	8.0	0.95	Yes	--
<i>Psophocarpus tetragonolobus</i> II	PTL II, WBA II	<i>Psophocarpus tetragonolobus</i> (Winged Bean) seeds	46	2	6.0	1.20	Yes	--
<i>Ricinus communis</i> I	RCA I, RCA ₁₂₀	<i>Ricinus communis</i> (Castor Bean) seeds	120	2	7.8	1.17	Yes	No
<i>Ricinus communis</i> II, ricin	RCA II, RCA ₆₀ , ricin	<i>Ricinus communis</i> (Castor Bean) seeds	60	1	7.1	1.17	Yes	No
Ricin A Chain	Ricin A Chain	RCA ₆₀	28	1	7.5	0.7	Yes	No
Ricin B Chain	Ricin B Chain	RCA ₆₀	32	1	4.5	1.64	Yes	No
<i>Sambucus nigra</i>	SNA, EBL	<i>Sambucus nigra</i> (Elderberry) bark	140	4	5.4 - 5.8	1.50	Yes	No
<i>Solanum tuberosum</i>	STL, PL	<i>Solanum tuberosum</i> , (potato) tubers	100	2	>9	0.80	Yes	No
<i>Sophora japonica</i>	SJA	<i>Sophora japonica</i> (Japanese Pagoda Tree) seeds	133	2	4.9 - 5.6	1.67	Yes	Ca ⁺⁺ , Mn ⁺⁺
Soybean	SBA	<i>Glycine max</i> (soybean) seeds	120	4	5.8 - 6	1.33	Yes	Ca ⁺⁺ , Mn ⁺⁺
<i>Ulex europaeus</i> I	UEA I	<i>Ulex europaeus</i> (Furze Gorse) seeds	63	2	4.5 - 5.1	1.30	Yes	Ca ⁺⁺ , Mn ⁺⁺ , Zn ⁺⁺
<i>Vicia villosa</i>	VVL, VVA	<i>Vicia villosa</i> (Hairy Vetch) seeds	144*	4	5.5 - 6.2	0.78	Yes	Ca ⁺⁺ , Mn ⁺⁺
Wheat Germ	WGA	<i>Triticum vulgaris</i> (wheat germ)	36	2	>9	1.46	No	Ca ⁺⁺
Succinylated Wheat Germ	Succinylated WGA	<i>Triticum vulgaris</i> (wheat germ)	36	2	<3	1.46	No	Ca ⁺⁺
<i>Wisteria floribunda</i>	WFA, WFL	<i>Wisteria floribunda</i> (Japanese Wisteria) seeds	116	4	5.2 - 5.8	0.89	Yes	--

* Literature values reported: 102 kDa - 144 kDa

Lectin	Mitogenic Activity	Blood Group Specificity	Preferred Sugar Specificity	Inhibitor or Eluting Sugar	Special Application
<i>Narcissus pseudonarcissus</i>	No	Rabbit	α Man	Mec α Man	Detection of beginning of apoptosis of human cell lines
Peanut	No	T antigen (M, N)	Gal β 3GalNAc	Gal	Detection of relocalization of Tag antigen in large bowel carcinoma
<i>Phaseolus vulgaris</i> Erythroagglutinin (PHA-E)	Yes	A(-SA)	Gal β 4GlcNAc β 2Man α 6 (GlcNAc β 4) (GlcNAc β 4Man α 3) Man β 4	bovine thyroglobulin, acetic acid	Binding to central nervous system cells (HNK-1 antigen)
<i>Phaseolus vulgaris</i> Agglutinin (PHA-E+L)	Yes	A(-SA)	See PHA-E/PHA-L	bovine thyroglobulin, acetic acid	Binding to central nervous system cells (HNK-1 antigen)
<i>Phaseolus vulgaris</i> Leucoagglutinin (PHA-L)	Yes	–	Gal β 4GlcNAc β 6(GlcNAc β 2Man α 3)Man α 3	bovine thyroglobulin, acetic acid	Anterograde neuronal tracing; metastatic tumor marker; lymphocyte mitogen for lymphokine production
<i>Pisum sativum</i>	Yes	Non-specific	α Man, α Glc	Mec α Man+ Mec α Glc	Separation of lymphoblastic leukemia antigen in kidney cells; purification of feline T-lymphocytes from peripheral blood
<i>Psophocarpus tetragonolobus</i> I	No	Rabbit, O(-SA)	GalNAc, Gal	GalNAc	Staining of mouse M-cells; blood vessel staining of A and B blood group individuals
<i>Psophocarpus tetragonolobus</i> II	No	O(-SA)	GalNAc, Gal	GalNAc	Blood vessel staining (human) of blood group O individuals
<i>Ricinus communis</i> I	No	Non-specific	Gal	Gal or Lactose	Labeling of receptors on sprouting rat neurons
<i>Ricinus communis</i> II, ricin	No	Non-specific	Gal, GalNAc	Gal or Lactose	TOXIC - Used in rat neuronal retrograde transport (suicide transport)
Ricin A Chain	No	–	–	–	Produce hybrid toxins
Ricin B Chain	Yes	–	Gal	Gal or Lactose	Used to potentiate antibody ricin A chain conjugates for tumor toxicity
<i>Sambucus nigra</i>	No	Non-specific	Neu5Ac α 6Gal/GalNAc	Lactose in buffered saline & acetic acid	Used to distinguish sialylated oligosaccharides bound by human A influenza virus
<i>Solanum tuberosum</i>	No	Non-specific	(GlcNAc) ₂₋₄	Chitin hydrolysate	Staining of prostate cancer cell line; bacterial cell wall binding
<i>Sophora japonica</i>	No	A>B>O(-SA)	β GalNAc	GalNAc	Distinguishes between pathogenic and non pathogenic trypanosomes
Soybean	Yes	A>O>B	α > β GalNAc	GalNAc	Stem cell separation
<i>Ulex europaeus</i> I	No	O>A2	α Fuc	L-Fuc	Human endothelial cell marker
<i>Vicia villosa</i>	No	Tn antigen	GalNAc	GalNAc	Staining of neurons on human cerebral cortex, detection of Tn and Cad antigens
Wheat Germ	Yes	A,B,O	GlcNAc	Chitin hydrolysate or GlcNAc with acid or salt	Insulin receptor purification, neuronal tracing; bacterial cell wall binding
Succinylated Wheat Germ	No	A,B,O	GlcNAc	Chitin hydrolysate or GlcNAc with acid or salt	Differential binding to intrahepatic blood vessels
<i>Wisteria floribunda</i>	Yes	Non-specific	GalNAc	GalNAc, acetic acid	Serotyping α -hemolytic streptococci

Sugar Abbreviations:

Fuc	L-Fucose	Man	Mannose
Gal	D-Galactose	Mec α Glc	α -Methylglucoside
GalNAc	N-Acetylgalactosamine	Mec α Man	α -Methylmannoside
Glc	D-Glucose	Neu5Ac	N-Acetylneuraminic acid (sialic acid)
GlcNAc	N-Acetylglucosamine	SA	Sialic Acid

Lectin Products

Agaricus Bisporus Lectin (ABL)

Unconjugated	L-1420	2 mg
Fluorescein	FL-1421	1 mg
Agarose (2 mg lectin/ml gel)	AL-1423	2 ml
Biotin	B-1425	1 mg

Aleuria Aurantia Lectin (AAL)

Unconjugated	L-1390	2 mg
Alkaline Phosphatase	MB-4100	150 µg
Agarose (2 mg lectin/ml gel)	AL-1393	2 ml
Biotin	B-1395	1 mg

Amaranthus Caudatus Lectin (ACL, ACA)

Unconjugated	L-1250	5 mg
Fluorescein	FL-1251	2 mg
Biotin	B-1255	2 mg

Bauhinia Purpurea Lectin (BPL)

Unconjugated	L-1280	5 mg
Fluorescein	FL-1281	2 mg
Biotin	B-1285	2 mg

Concanavalin A (Con A)

Unconjugated	L-1000	500 mg
Fluorescein	FL-1001	25 mg
Rhodamine	RL-1002	25 mg
Agarose (6 mg lectin/ml gel)	AL-1003	10 ml, 100 ml
Biotin	B-1005	5 mg

Succinylated Concanavalin A

Unconjugated	L-1000S	25 mg
Fluorescein	FL-1001S	10 mg
Agarose (6 mg lectin/ml gel)	AL-1003S	2 ml
Biotin	B-1005S	5 mg

Datura Stramonium Lectin

Unconjugated	L-1180	5 mg
Fluorescein	FL-1181	2 mg
Agarose (3 mg lectin/ml gel)	AL-1183	2 ml
Biotin	B-1185	2 mg

Dolichos Biflorus Agglutinin (DBA)

Unconjugated	L-1030	5 mg, 10 mg
Fluorescein	FL-1031	2 mg, 5 mg
Rhodamine	RL-1032	2 mg
Agarose (3 mg lectin/ml gel)	AL-1033	2 ml
Biotin	B-1035	5 mg

Erythrina Cristagalli Lectin (ECL, ECA)

Unconjugated	L-1140	10 mg
Fluorescein	FL-1141	5 mg
Agarose (3 mg lectin/ml gel)	AL-1143	2 ml
Biotin	B-1145	5 mg

Euonymus Europaeus Lectin (EEL)

Unconjugated	L-1330	5 mg
Fluorescein	FL-1331	2 mg
Biotin	B-1335	2 mg

Galanthus Nivalis Lectin (GNL)

Unconjugated	L-1240	5 mg
Fluorescein	FL-1241	2 mg
Agarose (3 mg lectin/ml gel)	AL-1243	2 ml, 5 ml
Biotin	B-1245	2 mg

Griffonia (Bandeiraea) Simplicifolia Lectin I (GSL I, BSL I)

Unconjugated	L-1100	5 mg
Fluorescein	FL-1101	2 mg, 5 mg
Rhodamine	RL-1102	2 mg
Agarose (4 mg lectin/ml gel)	AL-1103	2 ml
Biotin	B-1105	2 mg

GSL I – isolectin B₄

Unconjugated	L-1104	1 mg
Fluorescein	FL-1201	0.5 mg
DyLight® 594	DL-1207	0.5 mg
Biotin	B-1205	0.5 mg

Griffonia (Bandeiraea) Simplicifolia Lectin II (GSL II, BSL II)

Unconjugated	L-1210	5 mg
Fluorescein	FL-1211	2 mg
Agarose (3 mg lectin/ml gel)	AL-1213	2 ml
Biotin	B-1215	2 mg

Hippeastrum Hybrid (Amaryllis) Lectin (HHL, AL)

Unconjugated	L-1380	5 mg
Biotin	B-1385	2 mg

Jacalin

Unconjugated	L-1150	10 mg, 25 mg
Fluorescein	FL-1151	5 mg
Agarose (4 mg lectin/ml gel)	AL-1153	2 ml, 5 ml, 10 ml
Biotin	B-1155	5 mg

Lens Culinaris Agglutinin (LCA, LcH)

Unconjugated	L-1040	10 mg, 25 mg
Fluorescein	FL-1041	5 mg
Rhodamine	RL-1042	5 mg
Agarose (3 mg lectin/ml gel)	AL-1043	2 ml, 5 ml, 10 ml
Biotin	B-1045	5 mg

Lotus Tetragonolobus Lectin (LTL)

Unconjugated	L-1320	5 mg
Fluorescein	FL-1321	2 mg
Agarose (3 mg lectin/ml gel)	AL-1323	2 ml
Biotin	B-1325	2 mg

Lycopersicon Esculentum (Tomato) Lectin (LEL, TL)

Unconjugated	L-1170	2 mg
Fluorescein	FL-1171	1 mg
Texas Red®	TL-1176	1 mg
DyLight® 488	DL-1174	1 mg
DyLight® 594	DL-1177	1 mg
DyLight® 649	DL-1178	1 mg
Agarose (2 mg lectin/ml gel)	AL-1173	2 ml
Biotin	B-1175	1 mg

Maackia Amurensis Lectin I (MAL I, MAL)

Unconjugated	L-1310	5 mg
Fluorescein	FL-1311	2 mg
Biotin	B-1315	2 mg

Maackia Amurensis Lectin II (MAL II, MAH)

Unconjugated	L-1260	2 mg
Biotin	B-1265	1 mg

Maclura Pomifera Lectin (MPL)

Unconjugated	L-1340	5 mg
Fluorescein	FL-1341	2 mg
Biotin	B-1345	2 mg

Musa Paradisiaca Lectin (BanLec)

Unconjugated	L-1410	5 mg
Fluorescein	FL-1411	2 mg
Agarose (5 mg lectin/ml gel)	AL-1413	2 ml
Biotin	B-1415	2 mg

Narcissus Pseudonarcissus (Daffodil) Lectin (NPL, NPA, DL)

Unconjugated	L-1370	5 mg
Biotin	B-1375	2 mg

Peanut Agglutinin (PNA)

Unconjugated	L-1070	5 mg, 10 mg, 25 mg
Fluorescein	FL-1071	5 mg, 10 mg
Rhodamine	RL-1072	5 mg
Agarose (5 mg lectin/ml gel)	AL-1073	2 ml, 5 ml
Biotin	B-1075	5 mg

Phaseolus Vulgaris Agglutinin (PHA)

Phaseolus vulgaris Erythroagglutinin (PHA-E)

Unconjugated	L-1120	5 mg
Fluorescein	FL-1121	2 mg
Agarose (3 mg lectin/ml gel)	AL-1123	2 ml
Biotin	B-1125	2 mg

Phaseolus vulgaris Agglutinin (PHA-E+L)

Unconjugated	L-1220	25 mg
--------------	--------	-------

Phaseolus vulgaris Leucoagglutinin (PHA-L)

Unconjugated	L-1110	5 mg
Fluorescein	FL-1111	2 mg
Rhodamine	RL-1112	2 mg
Agarose (3 mg lectin/ml gel)	AL-1113	2 ml
Biotin	B-1115	2 mg

Pisum Sativum Agglutinin (PSA)

Unconjugated	L-1050	10 mg
Fluorescein	FL-1051	5 mg
Rhodamine	RL-1052	5 mg
Agarose (3 mg lectin/ml gel)	AL-1053	2 ml
Biotin	B-1055	5 mg

Psophocarpus Tetragonolobus Lectin I (PTL I, WBA I)

Unconjugated	L-1360	5 mg
Biotin	B-1365	2 mg

Lectin Products (continued)

Psophocarpus Tetragonolobus Lectin II (PTL II, WBA II)

Unconjugated	L-1400	5 mg
Biotin	B-1405	2 mg

Ricinus Communis Agglutinin I (RCA I, RCA₁₂₀)

Unconjugated	L-1080	10 mg
Fluorescein	FL-1081	5 mg
DyLight® 594	DL-1087	1 mg
Rhodamine	RL-1082	5 mg
Agarose (4 mg lectin/ml gel)	AL-1083	2 ml, 5 ml, 10 ml
Biotin	B-1085	5 mg

Ricinus Communis Agglutinin II (RCA II, RCA₆₀, ricin)

Unconjugated	L-1090	10 mg
Fluorescein	FL-1091	5 mg
Agarose (4 mg lectin/ml gel)	AL-1093	2 ml
Biotin	B-1095	5 mg

Ricin A Chain

Unconjugated	L-1190	1 mg
--------------	--------	------

Ricin B Chain

Unconjugated	L-1290	1 mg
--------------	--------	------

Sambucus Nigra Lectin (SNA, EBL)

Unconjugated	L-1300	5 mg
Fluorescein	FL-1301	2 mg
Agarose (3 mg lectin/ml gel)	AL-1303	2 ml
Biotin	B-1305	2 mg

Solanum Tuberosum (Potato) Lectin (STL, PL)

Unconjugated	L-1160	5 mg
Fluorescein	FL-1161	2 mg
Biotin	B-1165	2 mg

Sophora Japonica Agglutinin (SJA)

Unconjugated	L-1130	5 mg
Biotin	B-1135	2 mg

Soybean Agglutinin (SBA)

Unconjugated	L-1010	10 mg, 25 mg
Fluorescein	FL-1011	2 mg
Rhodamine	RL-1012	2 mg
Agarose (4 mg lectin/ml gel)	AL-1013	2 ml
Biotin	B-1015	5 mg

Ulex Europaeus Agglutinin I (UEA I)

Unconjugated	L-1060	2 mg, 5 mg
Fluorescein	FL-1061	2 mg, 5 mg
Rhodamine	RL-1062	2 mg
Agarose (2 mg lectin/ml gel)	AL-1063	2 ml
Biotin	B-1065	2 mg

Vicia Villosa Lectin (VVL, VVA)

Unconjugated	L-1230	5 mg
Fluorescein	FL-1231	2 mg
Agarose (3 mg lectin/ml gel)	AL-1233	2 ml
Biotin	B-1235	2 mg

Wheat Germ Agglutinin (WGA)

Unconjugated	L-1020	10 mg, 25 mg
Fluorescein	FL-1021	5 mg, 10 mg
Rhodamine	RL-1022	5 mg, 10 mg
Peroxidase	PL-1026	2 mg
Agarose (7 mg lectin/ml gel)	AL-1023	2 ml, 5 ml, 10 ml
Biotin	B-1025	5 mg

Succinylated Wheat Germ Agglutinin

Unconjugated	L-1020S	10 mg
Fluorescein	FL-1021S	5 mg
Rhodamine	RL-1022S	5 mg
Agarose (3 mg lectin/ml gel)	AL-1023S	2 ml, 5 ml
Biotin	B-1025S	5 mg

Wisteria Floribunda Lectin (WFA, WFL)

Unconjugated	L-1350	5 mg
Fluorescein	FL-1351	2 mg
Agarose (3 mg lectin/ml gel)	AL-1353	2 ml
Biotin	B-1355	2 mg

Lectin Screening Kits

Lectin Screening Kits I

Unconjugated Lectin Kit I	LK-2000	1 kit
Biotinylated Lectin Kit I	BK-1000	1 kit
Fluorescein Lectin Kit I	FLK-2100	1 kit
Rhodamine Lectin Kit I	RLK-2200	1 kit

Lectin Screening Kits II

Lectin Kit II	LK-3000	1 kit
Biotinylated Lectin Kit II	BK-2000	1 kit
Fluorescein Lectin Kit II	FLK-3100	1 kit
Rhodamine Lectin Kit II	RLK-3200	1 kit

Lectin Screening Kits III

Biotinylated Lectin Kit III	BK-3000	1 kit
Fluorescein Lectin Kit III	FLK-4100	1 kit

Kit I (LK-2000, BK-1000, FLK-2100, RLK-2200) consists of 1 mg of the following lectins or lectin conjugates: Con A, DBA, PNA, RCA I, SBA, UEA I, WGA.

Kit II (LK-3000, BK-2000, FLK-3100, RLK-3200) consists of 1 mg of the following lectins or lectin conjugates: GSL I, LCA, PHA-E, PHA-L, PSA, Succinylated WGA, SJA*.

*SJA is not included in FLK-3100

Kit III (BK-3000, FLK-4100) consists of 0.5 mg of the following lectin conjugates: DSL, ECL, GSL II, Jacalin, LEL, STL, VVL.

Inhibiting Sugars

Product	Catalog Number	Unit Size	Stock Concentration*
Chitin Hydrolysate	SP-0090	10 ml	N.A.
Sugars			
<i>N</i> -acetylgalactosamine	S-9001	111 mg	100 mM
<i>N</i> -acetylglucosamine	S-9002	442 mg	400 mM
galactose	S-9003	360 mg	400 mM
lactose	S-9004	721 mg	400 mM
α -methylmannoside	S-9005	388 mg	400 mM
α -methylglucoside	S-9006	388 mg	400 mM
L-fucose	S-9007	82 mg	100 mM
<i>N</i> -acetylneuraminic acid (sialic acid)	S-9008	619 mg	400 mM

* Stock concentration if reconstituted in 5 ml.

Antibodies to Lectins

Product	Conjugate	Catalog Number	Unit Size
Anti-Concanavalin A	Unconjugated	AS-2004	1 mg
Anti- <i>Dolichos biflorus</i> agglutinin	Unconjugated	AS-2034	1 mg
Anti- <i>Galanthus nivalis</i> lectin	Unconjugated	AS-2240	1 mg
Anti- <i>Griffonia (Bandeiraea) simplicifolia</i> lectin I	Unconjugated	AS-2104	1 mg
Anti- <i>Lens culinaris</i> agglutinin/ <i>Pisum sativum</i> agglutinin	Unconjugated	AS-2044	1 mg
Anti-Peanut agglutinin	Unconjugated	AS-2074	1 mg
	Biotinylated	BA-0074	0.5 mg
Anti- <i>Phaseolus vulgaris</i> agglutinin (E+L)	Unconjugated	AS-2224	1 mg
	Biotinylated	BA-0224	0.5 mg
Anti- <i>Phaseolus vulgaris</i> agglutinin (E+L)*	Unconjugated	AS-2300	1 mg
Anti- <i>Ricinus communis</i> agglutinin I & II	Unconjugated	AS-2084	1 mg
	Biotinylated	BA-0084	0.5 mg
Anti-Soybean agglutinin	Unconjugated	AS-2014	1 mg
Anti- <i>Ulex europaeus</i> agglutinin I	Unconjugated	AS-2064	1 mg
	Biotinylated	BA-0064	0.5 mg
Anti-Wheat Germ agglutinin	Unconjugated	AS-2024	1 mg
	Biotinylated	BA-0024	0.5 mg

Glycobiology Reagents

Neoglycoprotein

BSA-Mannose	G-1000	10 mg
BSA-Galactose	G-2000	10 mg
BSA-Fucose	G-3000	10 mg
BSA-GalNAc	G-4000	10 mg
BSA-Chitin	G-5000	10 mg
BSA-GlcNAc	G-6000	10 mg
BSA-Sialic Acid	G-7000	10 mg
BSA-Arabinose	G-8000	10 mg
BSA-Glucose	G-9000	10 mg

Neoglycoproteins generally describe proteins that, in their native form, do not contain carbohydrate groups but have been chemically derivatized with sugars.

Nine neoglycoproteins are available, each with a unique sugar attached to the protein. They are produced by coupling simple sugars to highly purified Bovine Serum Albumin (BSA). These glycoconjugates retain free amino groups, allowing them to be subsequently labeled with fluorochromes, biotin, or other tags. They can be used: 1) to check the binding activity of agarose lectin columns; 2) as a positive control for glycoprotein blots; 3) to localize endogenous lectins in tissues and cells. Each is supplied as a 10 mg/ml solution.

Agarose-Bound Lectins	Sugars Conjugated to Bovine Serum Albumin								
	Mannose (G-1000)	Galactose (G-2000)	Fucose (G-3000)	GalNAc (G-4000)	[GlcNAc]1-3 (G-5000)*	GlcNAc (G-6000)	Sialic Acid (G-7000)	Arabinose (G-8000)	Glucose (G-9000)
ABL (AL-1423)				+					
AAL (AL-1393)			+					+	
BanLec (AL-1413)	+								+
Con A (AL-1003)	+								+
sCon A (AL-1003S)	+								+
DSL (AL-1183)					+				
DBA (AL-1033)				+					
ECA (AL-1143)				+					
GNA (AL-1243)	+								
GSL I (AL-1103)		+							
GSL II (AL-1213)						+			
Jacalin (AL-1153)		+							
LCA (AL-1043)	+								+
LTL (AL-1323)			+					+	
LEL (AL-1173)					+				
PNA (AL-1073)		+							
PHA-E (AL-1123)		+							
PHA-L (AL-1113)		+							
PSA (AL-1053)	+								+
RCA₁₂₀ (AL-1083)		+							
RCA₆₀ (AL-1093)		+							
SNA (AL-1303)							+		
SBA (AL-1013)				+					
UEA I (AL-1063)			+					+	
VVA (AL-1233)				+					
WGA (AL-1023)						+			
sWGA (AL-1023S)						+			
WFL (AL-1353)				+					

† Indicates binding to lectin.

* [GlcNAc]₁₋₃ is predominant sugar mixture in chitin hydrolysate.

Glycoprotein Eluting Solution for Agarose Bound:

Mannose- or Glucose-binding Lectins	ES-1100	100 ml
Galactose- or GalNAc-binding Lectins	ES-2100	100 ml
Fucose- or Arabinose-binding Lectins	ES-3100	100 ml
GlcNAc- or Chitin-binding Lectins	ES-5100	100 ml
Sialic Acid-binding Lectins	ES-7100	100 ml

Glycoproteins are frequently isolated and purified from protein mixtures using columns of agarose-bound lectins. After applying a protein mixture, the agarose-lectin column is washed free of unwanted proteins and the glycoprotein bound to the lectin is eluted with a sugar that inhibits binding. Unfortunately, achieving complete elution with a simple sugar solution can be difficult. Vector Laboratories has developed five Glycoprotein Elution Solutions in the neutral pH range that maximize the yield of eluted glycoproteins and preserve the activity of the agarose-bound lectins for repeated use. All components of these ready-to-use Glycoprotein Eluting Solutions can subsequently be removed by dialysis.

Agarose-Bound Lectins	Glycoprotein Eluting Solutions				
	ES-1100	ES-2100	ES-3100	ES-5100	ES-7100
ABL (AL-1423) *					
AAL (AL-1393)			+		
BanLec (AL-1413)	+				
Con A (AL-1003)	+				
sCon A (AL-1003S)	+				
DSL (AL-1183)				+	
DBA (AL-1033)		+			
ECA (AL-1143)		+			
GNA (AL-1243)	+				
GSL I (AL-1103)		+			
GSL II (AL-1213)				+	
Jacalin (AL-1153)		+			
LCA (AL-1043)	+				
LTL (AL-1323)			+		
LEL (AL-1173)				+	
PNA (AL-1073)		+			
PHA-E (AL-1123)		+			
PHA-L (AL-1113)		+			
PSA (AL-1053)	+				
RCA₁₂₀ (AL-1083)		+			
RCA₆₀ (AL-1093)		+			
SNA (AL-1303)					+
SBA (AL-1013)		+			
UEA I (AL-1063)			+		
VVA (AL-1233)		+			
WGA (AL-1023)				+	
sWGA (AL-1023S)				+	
WFL (AL-1353)		+			

† Indicates recommendation for eluting glycoproteins from agarose-bound lectins.

* Requires more stringent conditions for elution.



Notes:

ORDERING

(800) 227-6666

(650) 697-3600

Fax: (650) 697-0339

USA MAIN OFFICE:

Vector Laboratories, Inc.
30 Ingold Road
Burlingame, CA 94010

Ordering: (800) 227-6666
Technical Service: (650) 697-3600

8:30 am to 5:00 pm Pacific Time

Fax: (650) 697-0339

Website: www.vectorlabs.com

Email: (USA) vector@vectorlabs.com
(Outside USA) intl@vectorlabs.com

ORDERING INFORMATION:

Orders may be placed by mail, telephone, fax, email or secure online ordering. Please include the following with each order:

- Product name and catalog number
- Unit size and quantity
- Billing and shipping addresses
- Purchase order number
- Name and phone number of person placing order

For the convenience of our customers, orders using VISA, Mastercard, or American Express are accepted and processed immediately. Telephone orders over \$2000 may require written/faxed confirmation. A confirmation should be boldly marked "Confirming Order. Do Not Duplicate." Duplicate shipments due to incorrectly marked confirming orders cannot be returned for credit. **No returned product will be accepted or credited without prior authorization from Vector Laboratories.**

DISCOUNT POLICY:

Discounts for most products are available for large quantities or bulk orders. Orders for 10-29 of the same product size will automatically receive a 10% discount, orders for 30 or more will receive a 20% discount. Please inquire for pricing of bulk quantities of any product. All RCA products are excluded from our regular discount schedule.

PAYMENT/SHIPPING TERMS:

Our payment terms are net 30 days. Prices are FCA Burlingame, California. Shipping charges will be prepaid and added to the invoice. Orders are usually shipped the same day they are received. Unless requested otherwise, all products are shipped 2nd Day Air. RCA₆₀ products are shipped in the USA only according to federal transportation regulations requiring additional shipping charges. We require a written confirmation for all RCA₆₀ products.

INTERNATIONAL ORDERS:

International orders are shipped by air. For the convenience of our customers in the U.K. and Canada, Vector Laboratories maintains a complete inventory of all Vector products at our subsidiaries:

United Kingdom Office

Vector Laboratories Ltd.,
3 Accent Park, Bakewell Road, Orton Southgate,
Peterborough, PE2 6XS, United Kingdom
tel (01733) 237999, fax (01733) 237119
email: vector@vectorlabs.co.uk

Canadian Office

Vector Laboratories (Canada), Inc.
3390 South Service Road
Burlington, Ontario, L7N 3J5, Canada
tel (905) 681-0900, fax (905) 681-0166
tel (888) 629-2121 (orders only)
email: vectorca@netinc.ca

Vector products are also available from over 45 distributors worldwide. Please contact us or visit our website for a complete listing of these representatives.

Although provided in a highly purified form, our products are not intended for clinical diagnosis or drug use, nor have they been packaged under sterile conditions. All products listed in this catalog are for research purposes only. The listing of any product in this catalog does not imply the absence of a patent covering its use, does not constitute license under any existing or pending patent, nor is it intended or implied as a recommendation for the use of such products in infringement of any patent. The responsibility for determining the existence of such patents rests solely with the user.

Trademarks of Vector: ABC-Amp, Animal-Free Blocker, AT-Yellow, BLOXALL, DuoLuX, Elite, EndTag, FastTag, Fusion Aid, GC-Red, HYBEX, ImmEdge, ImmPACT, ImmPRESS, ImmPrint, M.O.M., Neurobiotin, NickKit p.s.o., Photoprobe, PolyBlock, ProFusion, ProtOn, Quant*Tag, Resolve-It, UltraSNAP, Vectabond, VectaCell, VectaMount, Vectashield, Vectashield Hard*Set, Vectastain, Vector, Vector Black, Vector Blue, Vector NovaRED, Vector Red, Vectrex. Others: ABTS-Boehringer Mannheim; Texas Red-Molecular Probes; DyLight-Thermo Fisher Scientific; APS-5-Lumigen Inc. DuoLuX™ is a special formula of APS-5 produced for Vector Labs by Lumigen Inc. Southfield, MI.



This catalog has been printed on recycled paper.



Printed in the USA

©Vector Laboratories, Inc. 2016





VECTOR LABORATORIES INC.

30 Ingold Road, Burlingame, CA 94010 u.s.a.

TEL: (650) 697-3600 FAX: (650) 697-0339

www.vectorlabs.com