
From Eye to Insight








In molecular and genetic cell research, it's critical to witness transgene expression in early stages, to select the right sample to base your studies on. That's why you need a stereo microscope with superb brightness and a manual zooming option for fast and efficient screening. The Leica M205 FCA fluorescence stereo microscope now unites two worlds. It combines high-end imaging with super-fast screening to give you bright images fast.

Get a bright fluorescence signal


With separated, yet synchronized beam paths, one for the fluorescence excitation light and two for observation, the M205 FCA delivers an evenly and fully illuminated field of view at any zoom position you choose. The TripleBeam technology eliminates reflections in the observation path and leaves you a strong fluorescence signal with a noise-free, solid black background.

Finest details in 3D


FusionOptics technology overcomes optical limitations. Uniting high resolution and depth of field (DOF) it provides increased image brightness and facilitates orientation in 3D.



Solea senegalensis larvae nervous system, max. projection of a tile scan of 6 fields \times 33 planes. Parallax correction and tiling performed in LAS X after deconvolution with Huygens professional. Dr. Marco A. Campinho, CCMAR — Centre for Marine Sciences, Universidade do Algarve, Portugal.



With FusionOptics higher resolution can be achieved (red line) than with traditional optics (green line). At the same time more DOF is achieved.




2-in-1: the screening and imaging solution

With the large zoom range of 20.5:1 you can quickly change from a broad overview to the finest details. While the 1.0× PlanApo objective with a numerical aperture (NA) of 0.17 provides a large working distance the 2.0× CORR objective offers an unmatched NA of 0.35. Zoom-in manually and trust the fully coded system auto-storing all paramaters with your image for reliable and publishable results.

Continue working

Eliminate interruptions of your workflow: The unique four position coded filter changer allows you to work with several different fluorochromes. Move filters manually or automated and just continue working. Take advantage of the optional footswitches to easily change filters and focus or change the illumination setting during screening.







Leica MICROSYSTEMS

SPECIFICATIONS

LEICA M205 FCA

	EE137 (10/200) 07 (
OPTICAL DATA	20.5 : 1 manual, coded with FusionOptics
Zoom	
Data with standard optics (1× objective/10× eyepieces) - Zoom range - Resolution - Working distance - Object field	7.8×−160× max. 525 lp/mm 61.5 mm (planapochromatic) Ø 29.5 mm−1.44 mm
Maximum values (based on optics combination) - Magnification - Resolution - Visible structural width - Numerical aperture - Object field	1.280× 1,050 lp/mm 476 nm 0.35 Ø 59 mm
Working distances	135 mm (0.5× planachromatic) 112 mm (0.8× planachromatic) 67 mm (0.63× planapochromatic) 61.5 mm (1× planapochromatic) 30.5 mm (1.6× planapochromatic) 20.1 mm (2× planapochromatic)
OPTICS CARRIER	
100 % apochromatic optical system	CMO (Common Main Objective) lead-free
Specific surface resistivity (housing)	$2\times10^{\rm 11}~\Omega/mm^2$ discharge time <2 seconds from 1,000 V to 100 V
Encoded/automated	Zoom, iris diaphragm, objective nosepiece, filter
Engageable zoom notches	14 for repetitive tasks
Double-iris diaphragm for depth of field control	Built-in and encoded
ARTICLE NUMBER	10 450 826

Leica Microsystems (Schweiz) AG \cdot Max-Schmidheiny-Strasse 201 \cdot 9435 Heerbrugg, Switzerland T +41 71 726 34 34 \cdot F +41 71 726 34 44

www.leica-microsystems.com

CONNECT WITH US!

