



16:1
6:1

LEICA Z6 APOA

Leica

The image shows a close-up of a white Leica microscope objective lens. The lens is mounted on a black base. The text '16:1' and '6:1' is printed in large, bold, black font on the left side of the lens. The Leica logo is visible on the bottom of the lens. The text 'LEICA Z6 APOA' is printed vertically on the right side of the lens. The background is a plain, light-colored surface.

Leica Zoom Systems

Technical Information

Living up to Life

Leica
MICROSYSTEMS

Leica Z6 APO & Z16 APO, Z6 APO A & Z16 APO A

Technical Data, Performance Characteristics



Zoom system	Leica Z6 APO	Leica Z16 APO	Leica Z6 APO A	Leica Z16 APO A
Type	Apochromatic zoom system with central beam path, lead-free	Apochromatic zoom system with central beam path, lead-free	Apochromatic zoom system with central beam path, lead-free, motorized functions	Apochromatic zoom system with central beam path, lead-free, motorized functions
Zoom	6.3:1	16:1	6.3:1, motorized	16:1, motorized
Zoom factor	0.57× – 3.6×	0.57× – 9.2×	0.57× – 3.6×	0.57× – 9.2×
Built-in iris diaphragm	for adjusting the depth of field	for adjusting the depth of field	for adjusting the depth of field, motorized	for adjusting the depth of field, motorized
Switchable approachable positions for repetitive tasks	0.57 / 0.8 / 1 / 1.25 / 1.6 / 2 / 2.5 / 3.2 / 3.6	0.57 / 0.8 / 1 / 1.25 / 1.6 / 2 / 2.5 / 3.2 / 4 / 5 / 6.3 / 8 / 9.2	continuously motorized and approachable positions	continuously motorized and approachable positions
Zoom travel speed	–	–	1.6 s for zoom range	2.5 s for zoom range
Visual data with 1× planapochromatic objective / 10× eyepieces / 1.25× Y tube				
Magnification	7.1× – 45×	7.1 – 115×	7.1× – 45×	7.1× – 115×
Resolution	60 – 351 Lp/mm	51 – 336 Lp/mm	60 – 351 Lp/mm	51 – 336 Lp/mm
Visible structural width	8.3 – 1.4 μm	9.8 – 1.49 μm	8.3 – 1.4 μm	9.8 – 1.49 μm
Numerical aperture	0.02 – 0.117nA	0.017 – 0.112nA	0.02 – 0.117nA	0.017 – 0.112nA
Object field Ø	32.3 mm – 5.1 mm	32.3 mm – 2.0 mm	29.5 mm – 4.7 mm	29.5 mm – 1.83 mm
Depth of field (diaphragm open)	3.1 mm – 0.09 mm	3.8 mm – 0.05 mm	3.1 mm – 0.09 mm	3.8 mm – 0.05 mm
Depth of field (diaphragm closed)	18.1 mm – 0.4 mm	18.4 mm – 0.4 mm	18.1 mm – 0.4 mm	18.4 mm – 0.4 mm
Visual data with 2× planapochromatic objective / 40× eyepieces / 1.25× Y tube				
Magnification	57× – 360×	57× – 920×	57× – 360×	57× – 920×
Resolution	120 – 702 Lp/mm	102 – 672 Lp/mm	120 – 702Lp/mm	102 – 672Lp/mm
Visible structural width	4.2 – 0.7 μm	4.9 – 0.74 μm	4.2 – 0.7 μm	4.9 – 0.74 μm
Numerical aperture	0.04 – 0.234nA	0.034 – 0.224nA	0.04 – 0.234nA	0.034 – 0.224nA
Object field Ø	4.2 mm – 0.67 mm	4.2 mm – 0.26 mm	4.2 mm – 0.67 mm	4.2 mm – 0.26 mm
Data with Leica DFC490 digital camera / 1× planapochromatic objective / AS tube / 0.63× video objective				
Magnification Chip: Object	0.36× – 2.3×	0.36× – 5.8×	0.36× – 2.3×	0.36× – 5.8×
Digital resolution*	33.3 – 210 Lp/mm	33.3 – 336 Lp/mm	33.3 – 210Lp/mm	33.3 – 336Lp/mm
Object field projected onto chip	24.5 mm × 18.4 mm / 3.9 mm × 2.9 mm	24.5 mm × 18.4 mm / 1.5 mm × 1.14 mm	24.5 mm × 18.4 mm / 3.9 mm × 2.9 mm	24.5 mm × 18.4 mm / 1.5 mm × 1.14 mm
Depth of field (diaphragm open)	1.06 mm – 0.03 mm	1.4 mm – 0.03 mm	1.06 mm – 0.03 mm	1.4 mm – 0.03 mm
Depth of field (diaphragm closed)	10.7 mm – 0.26 mm	10.9 mm – 0.3 mm	10.7 mm – 0.26 mm	10.9 mm – 0.3 mm



Optical accessories	Leica Z6 APO & Z6 APO	Leica Z6 APO A & Z16 APO A
Objectives	<ul style="list-style-type: none"> – planapochromatic 1×, 2×, 0.8×, 0.5×, 5× – achromatic objectives M series 0.63×, 0.5×, 0.32×, lead-free 	<ul style="list-style-type: none"> – planapochromatic 1×, 2×, 0.8×, 0.5×, 5× – achromatic objectives M series 0.63×, 0.5×, 0.32×, lead-free
Working distance	<ul style="list-style-type: none"> – 187 mm (planapochromat 0.5×) – 97 mm (planapochromat 1×) – 112 mm (planapochromat 0.8×) – 39 mm (planapochromat 2×) – 19 mm (planapochromat 5.0×) 	<ul style="list-style-type: none"> – 187 mm (planapochromat 0.5×) – 97 mm (planapochromat 1×) – 112 mm (planapochromat 0.8×) – 39 mm (planapochromat 2×) – 19 mm (planapochromat 5.0×)
Objective adapter	<ul style="list-style-type: none"> – for M series achromats – for HR objectives 10× and 20× – for DM objectives 10× and 20× 	<ul style="list-style-type: none"> – for M series achromats – for HR objectives 10× and 20× – for DM objectives 10× and 20×
HR objectives*	<ul style="list-style-type: none"> – HR 10×/0.45, working distance 19 mm – HR 20×/0.42, working distance 13 mm 	<ul style="list-style-type: none"> – HR 10×/0.45, working distance 19 mm – HR 20×/0.42, working distance 13 mm
DM objectives*	<ul style="list-style-type: none"> – DM objective N Plan L 20× /0.40 corr – DM objective N Plan 10× /0.25–/A5.8 	<ul style="list-style-type: none"> – DM objective N Plan L 20× /0.40 corr – DM objective N Plan 10× /0.25–/A5.8
Fine focusing	10-mm path, optional	10-mm path, motorized and integrated
Binocular tubes, ergonomics	<ul style="list-style-type: none"> – inclined and straight binocular tubes – apochromatic ErgoTube® 10° – 50° with synchronous interpupillary distance adjustment – various ErgoModules® (optional) 	<ul style="list-style-type: none"> – inclined and straight binocular tubes – apochromatic ErgoTube® 10° – 50° with synchronous interpupillary distance adjustment – various ErgoModules® (optional)
	<i>ErgoTube® and ErgoModule® are registered in the United States Patent and Trademark Office</i>	<i>ErgoTube® and ErgoModule® are registered in the United States Patent and Trademark Office</i>
Interpupillary distance	55 mm – 75 mm	55 mm – 75 mm
Ergonomic wide-field eyepieces for persons wearing glasses	10×/21, 16×/14, 25×/9.5, 40×/6, distortion-free plug-on eyecups to protect against infections	10×/21, 16×/14, 25×/9.5, 40×/6, distortion-free plug-on eyecups to protect against infections
Electrical Interface		RS232, USB via motorized focus

* suitable for upper zoom range only



Optical Data – Visual Data with Y Tube 1.25×

Objectives		1× planapochromat		0.5× planapochromat		0.8× planapochromat		2× planapochromat		5× planapochromat	
Working distances		97 mm		187 mm		112 mm		39 mm		20 mm	
Eyepieces	Zoom position	Total mag.	Object field Ø	Total mag.	Object field Ø	Total mag.	Object field Ø	Total mag.	Object field Ø	Total mag.	Object field Ø
		×	mm	×	mm	×	mm	×	mm	×	mm
10×/23B	0.57	7.1	32.3	3.6	64.6	5.7	40.4	14.3	16.1	35.6	6.5
	0.8	10	23.0	5	46.0	8	28.8	20	11.5	50	4.6
	1	12.5	18.4	6.3	36.8	10	23.0	25	9.2	62.5	3.7
	1.25	15.6	14.7	7.8	29.4	12.5	18.4	31.3	7.4	78.1	2.9
	1.6	20	11.5	10	23.0	16	14.4	40	5.8	100	2.3
	2	25	9.2	12.5	18.4	20	11.5	50	4.6	125	1.8
	2.5	31.3	7.4	15.6	14.7	25	9.2	62.5	3.7	156.3	1.5
	3.2	40	5.8	20	11.5	32	7.2	80	2.9	200	1.2
	3.6	45	5.1	22.5	10.2	36	6.4	90	2.6	225	1.0
	4	50	4.6	25	9.2	40	5.8	100	2.3	250	0.9
	5	62.5	3.7	31.3	7.4	50	4.6	125	1.8	312.5	0.7
	6.3	78.8	2.9	39.4	5.8	63	3.7	157.5	1.5	393.8	0.6
8	100	2.3	50	4.6	80	2.9	200	1.2	500	0.5	
9.2	115	2.0	57.5	4.0	92	2.5	230	1.0	575	0.4	
16×/15B	0.57	11.4	21.1	5.7	42.1	9.1	26.3	22.8	10.5	57	4.2
	0.8	16	15.0	8	30.0	12.8	18.8	32	7.5	80	3.0
	1	20	12.0	10	24.0	16	15.0	40	6.0	100	2.4
	1.25	25	9.6	12.5	19.2	20	12.0	50	4.8	125	1.9
	1.6	32	7.5	16	15.0	25.6	9.4	64	3.8	160	1.5
	2	40	6.0	20	12.0	32	7.5	80	3.0	200	1.2
	2.5	50	4.8	25	9.6	40	6.0	100	2.4	250	1.0
	3.2	64	3.8	32	7.5	51.2	4.7	128	1.9	320	0.8
	3.6	72	3.3	36	6.7	57.6	4.2	144	1.7	360	0.7
	4	80	3.0	40	6.0	64	3.8	160	1.5	400	0.6
	5	100	2.4	50	4.8	80	3.0	200	1.2	500	0.5
	6.3	126	1.9	63	3.8	100.8	2.4	252	1.0	630	0.4
8	160	1.5	80	3.0	128	1.9	320	0.8	800	0.3	
9.2	184	1.3	92	2.6	147.2	1.6	368	0.7	920	0.3	
25×/9.5B	0.57	17.8	13.3	8.91	26.7	14.3	16.7	35.6	6.7	89	2.7
	0.8	25	9.5	12.5	19.0	20	11.9	50	4.8	125	1.9
	1	31.3	7.6	15.6	15.2	25	9.5	62.5	3.8	156	1.5
	1.25	39.1	6.1	19.5	12.2	31.3	7.6	78.1	3.0	195	1.2
	1.6	50	4.8	25	9.5	40	5.9	100	2.4	250	1.0
	2	62.5	3.8	31.3	7.6	50	4.8	125	1.9	313	0.8
	2.5	78.1	3.0	39.1	6.1	62.5	3.8	156	1.5	391	0.6
	3.2	100	2.4	50	4.8	80	3.0	200	1.2	500	0.5
	3.6	113	2.1	56.3	4.2	90	2.6	225	1.1	562	0.4
	4	125	1.9	62.5	3.8	100	2.4	250	1.0	625	0.4
	5	156	1.5	78.1	3.4	125	1.9	313	0.8	781	0.3
	6.3	197	1.2	98.4	2.4	158	1.5	394	0.6	984	0.2
8	250	1.0	125	1.9	200	1.2	500	0.5	1250	0.2	
9.2	288	0.8	144	1.6	230	1.0	575	0.4	1438	0.2	
40×/6B	0.57	28.5	8.4	14.3	16.8	22.8	10.5	57	4.2	143	1.7
	0.8	40	6.0	20	12.0	32	7.5	80	3.0	200	1.2
	1	50	4.8	25	9.6	40	6.0	100	2.4	250	1.0
	1.25	62.5	3.8	31.3	7.7	50	4.8	125	1.9	313	0.8
	1.6	80	3.0	40	6.0	64	3.8	160	1.5	400	0.6
	2	100	2.4	50	4.8	80	3.0	200	1.2	500	0.5
	2.5	125	1.9	62.5	3.8	100	2.4	250	1.0	625	0.4
	3.2	160	1.5	80	3.0	128	1.9	320	0.8	800	0.3
	3.6	180	1.3	90	2.7	144	1.7	360	0.7	900	0.3
	4	200	1.2	100	2.4	160	1.5	400	0.6	1000	0.2
	5	250	1.0	125	1.9	200	1.2	500	0.5	1250	0.2
	6.3	315	0.8	158	1.5	252	1.0	630	0.4	1575	0.2
8	400	0.6	200	1.2	320	0.8	800	0.3	2000	0.1	
9.2	460	0.5	230	1.0	368	0.7	920	0.3	2300	0.1	

Z6 APO / Z6 APO A: Zoom 0.57× – 3.6×
Z16 APO / Z16 APO A: Zoom 0.57× – 9.2×

Major User Segments

General	
Documentation	<ul style="list-style-type: none"> – Parallax-free precise imaging – Ideal for multi-focus applications
Polarization	– Axial imaging for true polarization colors
Measuring, analysis	– Parallax-free precise imaging
Technology, industry	
OEM	Checking the production flow, detection of position and orientation
Semiconductor industry	<ul style="list-style-type: none"> – Testing, checking, sorting semiconductor components – Inspecting gallium arsenide wafers (GaAs)
Microelectronics	Inspecting flex circuits and inkjet print heads
Fiber optic cable technology	<ul style="list-style-type: none"> – Aligning fiber optic cables – Inspecting laser diodes
Measuring and test facilities	Materials testing and measuring, placement of test probes
Metalworking industry	Quality control, testing metallographic specimens, documentation
Plastics industry	<ul style="list-style-type: none"> – Structural examinations of thin sections – Polarization-optical testing
Gemology	Inspection and classification of diamonds and gemstones
Implants	Control of coated stents (vascular supports) and catheters
Criminology	Sorting and documentation of trace materials such as fibers, human and animal hair, ammunition and documents (counterfeit). Documentation of impressions.
Natural science, life science	
Biology	Examination and documentation of illnesses, pests and ecological damages to plants
Geology	Testing petrological thin sections in polarized light
Medicine	Examination and documentation of histological thin sections
Education, training	Observing and demonstrating processes
Entomology	Documentation of beetles and insects
Pathology	Documenting organs



Name	MPixel	Sensor	Size	Pixel-Size	Cooling	Interface
DFC295	3	CMOS-c	1/2"	3.2µm	no	FWb
DFC290 HD	3	CMOS-c	1/2"	3.2µm	no	FWb+HDMI
DFC400	1.4	CCD-c	1/2"	4.65µm	no	FWb+Sync
DFC425	5	CCD-c	1/2"	2.78µm	no	FWb+Sync
DFC425 C	5	CCD-c	1/2"	2.78µm	yes	FWb+Sync
DFC495	8	CCD-c	2/3"	2.7µm	yes	FWb+Sync
DFC310 FX	1.4	CCD-c	2/3"	6.45µm	yes	FWb+Sync
DFC345 FX	2	CCD-m	1/2"	4.4µm	yes	FWb+Sync
DFC360 FX	1.4	CCD-m	2/3"	6.45µm	yes	FWb+Sync
DFC500	12	CCD-c	2/3"	6.45µm	yes	FWa+Sync

-c: color / -m: monochrome V.2010

Digital camera systems

- 12 730 209 Leica DFC295 camera kit
- 12 730 202 Leica DFC290 HD camera kit*
- 12 730 203 Leica DFC400 camera kit*
- 12 730 221 Leica DFC425 camera kit*
- 12 730 221 Leica DFC425 C camera kit*
- 12 730 223 Leica DFC495 camera kit*
- 11 547 001 Leica DFC310 FX camera kit for Mac
- 11 547 002 Leica DFC310 FX camera kit for PC
- 11 547 003 Leica DFC345 FX camera kit*
- 11 547 000 Leica DFC360 FX camera kit
- 12 730 066 Leica DFC500 camera kit
- 12 730 216 Leica IC80 HD camera (incl. USB cable and LAS EZ software)

- 12 730 228 Stand-alone kit (power pack, video cable, HDMI-cable, SD card, HD RC remote control)
- 12 730 229 Palm or foot switch with 2m cable

*no Mac-Support

Digital camera system accessories

- 12 446 120 FireWire cable, FWa, 3m, 6-pin to 6-pin
- 12 446 140 FireWire cable, FWa, 4m, 6-pin to 6-pin
- 12 730 186 FireWire cable, 3m, FWb, 9-pin to 9-pin
- 12 730 187 FireWire cable, FWab, 3m, 6-pin to 9-pin
- 12 730 211 HDMI cable, 3m, HDMI both ends
- 12 730 183 OHCI FireWire PCI32 card, FWa
- 12 730 225 PowerKit for DFC290 HD
- 12 730 210 OHCI FireWire PCI express card, FWb
- 12 447 066 Laptop PCMCIA FireWire interface card for FWa
- 12 730 188 Laptop power kit, for use with 4-pin or unpowered 6-pin FireWire ports

Photo tubes and C-mounts

- 10 447 319 Adapter T2, Canon EOS
- 10 447 436 DSLR-tube with T2 thread
- 10 446 175 SLR projection lens 2.5, for reflex cameras with T2 thread singlelens reflex cameras on video/photo tubes
- 10 445 928 0.32x video objective with C-mount for 1/3" digital cameras
- 10 445 929 0.5x video objective with C-mount for 1/2" digital cameras
- 10 447 367 0.63x video objective with C-mount for 2/3" digital cameras
- 10 446 307 0.8x video objective with C-mount for 2/3" digital cameras
- 10 445 930 1.0x video/photo objective
- 10 450 317 C-mount-adaptor for HC

ARTICLE DESCRIPTIONS

Eyepieces and graticules

Wide-field eyepieces for persons wearing glasses 10x /23, distortion-free, dioptic correction, with eyecups* – can lead to vignetting when used with other accessories, such as ergonomic or documentation tubes.

Remedy: Use ErgoTube 45° or 16x eyepieces

10 450 023

Wide-field eyepieces for persons wearing glasses 16x/15, distortion-free, dioptic correction, with eyecups*

Wide-field eyepieces for eyeglass wearers 25x/9.5, distortion-free, dioptic correction, with eyecups*

Wide-field eyepieces for eyeglass wearers 40x/6, distortion-free, dioptic correction, with eyecups*
* replaceable

Replacement eyecups for wide-field eyepieces for persons wearing glasses (not shown)

10 450 054

Graticule for encoded M series

10 450 116

Graticule, 5mm/0.05mm

10 450 117

Graticule, 10mm/0.1mm

10 450 118

Graticule, 100 Div./0.001

10 450 119

Graticule, crosshair

Tubes, ergonomics and integrated camera modules

Inclined binocular tube, 45° viewing angle

10 450 252

Binocular ErgoTube® with 45° viewing angle and long tubes, extension factor 1.6x

10 450 156

Straight binocular tube with 90° viewing angle

10 450 157

Binocular ErgoTube® with variable viewing angle 10° – 50°

10 450 042

Trinocular video-/phototube 50%, 30° viewing angle

10 450 043

Trinocular video-/phototube 100%, 30° viewing angle

10 450 044

Trinocular video-/phototube 100%, 5° – 45° viewing angle

10 446 123

ErgoWedge® 5° – 25° (not to be used with trinocular tubes)

10 346 910

ErgoWedge® ±15 (not to be used with trinocular tubes)

10 450 301

Video-/phototube HD F, 50%, 50%

10 450 302

Video-/phototube HD V, 100%, 50%, 50%, 100%

12 730 216

Leica IC80 HD camera kit

12 730 054

Leica IC D camera kit

10 447 158

Filter slide housing with 2 filter slides

Microscope carriers and accessories

10 446 343

Focusing arm for probers

10 450 174

Tiltable focusing arm

10 446 344

Focusing arm for columns Ø 25mm to bonders

10 447 259

Photo column, inclinable, Ø 25mm

10 450 173

Microscope carrier for focusing drive 10 450 217

10 450 106

Microscope carrier, compatible with M50/M80, S and Z series

10 447 109

Y tube for tubes of the Leica M series, tube factor 1.25x, light distribution 50/50

10 447 128

A tube

10 447 196


Carrier to AS tube at focusing drive

10 447 204

AS tube

10 447 170

Coaxial incident light housing for fiber-optic light



10 450 023

10 450 024

10 450 025

10 450 026

10 450 196

10 450 054

10 450 116

10 450 117

10 450 118

10 450 119

10 450 252

10 450 156

10 450 157

10 450 042

10 450 043

10 450 044

10 446 123

10 346 910

10 450 301

10 450 302

12 730 216

12 730 054

10 447 158

10 446 343

10 450 174

10 446 344

10 447 259

10 450 173

10 450 106

10 447 109

10 447 128

10 447 196

10 447 204

10 447 170

Zoom systems and accessories


- 10 447 174 Leica Z6 APO 6.3:1 apochromatic zoom system
- 10 447 173 Leica Z16 APO 16:1 apochromatic zoom system
- 10 446 368 Leica Z6 APO A 6.3:1 motorized zoom system, apochromatic, with distribution box and M12 cable
- 10 446 369 Leica Z16 APO A 16:1 motorized zoom system, apochromatic, with distribution box and M12 cable
- 10 447 458 Z16 APO A with objective slide
- 10 450 094 Objective slide with quarter-wave plate for Z6 APO / Z16 APO (please order objectives separately)
- 10 450 138 Objective slide for Z6 APO / Z16 APO (please order objectives separately)

Objectives and optical accessories



- 10 447 175 Fine focusing (for Z6 APO/Z16 APO only), 10 mm travel distance
- 10 447 177 0.5x planapochromatic objective
- 10 446 360 0.8x planapochromatic objective
- 10 447 176 1x planapochromatic objective
- 10 447 178 2x planapochromatic objective
- 10 447 243 Planapochromatic objective 5x
- 10 450 332 Adapter for M series achromatic objectives
- 10 450 191 0.32x achromatic objective
- 10 450 192 0.5x achromatic objective
- 10 450 160 0.63x achromatic objective
- 10 450 243 Adapter for LWD achromatic objectives
- 10 411 597 Achromatic objective f = 100 mm
- 10 441 787 Achromatic objective f = 150 mm
- 10 431 692 Achromatic objective f = 175 mm
- 10 382 162 Achromatic objective f = 200 mm
- 10 382 172 Achromatic objective f = 400 mm
- 35 000 011 Achromatic objective f = 800 mm
- 35 000 013 Achromatic objective f = 1200 mm
- 10 447 206 Adapter HR objectives to Z6/Z16
- 10 447 179 Adapter DM objectives to Z6/Z16
- 11 506 057 DM objective N plan L 20x/0.40 corr.*
- 11 506 084 DM objective N plan 10x/0.25 - /A5.8 *

*strong vignetting in the lower zoom range

- 10 450 169 LED3000 NV1™ – vertical illumination for M50/M80 for using objectives with major diameter 58 mm
- 10 450 266 Power supply unit for LED3000/LED5000
- 35 000 031 Protective glass for planapochromatic objectives
- 10 315 306 Analyzer in rotating mount
- 10 447 207 Quarter-wave plate, Ø 58 mm
- 10 447 460 Quarter-wave plate for 2.0x and 5.0x



4 Connection to Microscope carriers



Incident- and transmitted- light bases

- 10 447 342 Incident-light base medium
- 10 450 049 Incident-light base large
- 10 450 123 Transmitted-light base TL ST
- 10 450 124 Transmitted-light base BFDF
- 10 450 125 Transmitted-light base TL RC™ for external cold light sources
- 10 450 126 Transmitted-light base TL RC™ with integrated halogen illumination

Stages

- 10 447 269 Standard stage for TL BFDF, TL RC™ and TL RC™
- 10 450 127 Leica IsoPro™ manual mechanical stage for TL BFDF, TL RC™, TL RC™ transmitted-light bases and incident-light base (with adapter 10 450 122)
- 10 450 218 Leica IsoPro™ motorized mechanical stage for TL BFDF, TL RC™, TL RC™ transmitted-light bases and incident-light base (with adapter 10 450 122)
- 10 450 122 Adapter between mechanical stage and incident-light base

Leica MATS TL heating stage insert with control unit for TL transmitted-light bases


- 10 447 275 Adapter for stages with Ø120mm
- 10 447 391 Stage for LifeOnStage accessories
- 10 447 392 Universal carrier for Petri dishes, specimen slides (up to four) etc.

Gliding stage, Ø120 mm

- 10 446 301 Gliding stage, Ø120 mm
- 10 446 302 Polarization stage, Ø120mm
- 10 382 130 Attachable mechanical stage for polarization stage
- 10 361 719 Sensitive-tint plate for Pol rotating stage
- 10 446 303 Cup stage, Ø120 mm
- 10 446 228 Glass stage plate with Pol, Ø120 mm
- 10 450 058 Stage plate, b/w for TL bases
- 10 450 059 Replacement knobs for IsoPro™ manual mechanical stage

Focusing drives

- 10 450 172 Coarse focusing drive with profile column 500mm
- 10 450 257 Focusing drive, coarse/fine, with 620mm profile column
- 10 450 258 Motorized focus with profile column 620 mm
- 10 450 289 XL extension – for viewing large specimens
- 10 450 265 Arm for fastening LED spotlights
- 10 450 205 Combination light guide adapter for installation on the focusing column
- 10 447 369 CTL2 cable, 2 m



Filters

- 10 447 400 Daylight filter for TL ST base
- 10 447 394 BG38 filter for TL RC™/RCI™ transmitted-light base
- 10 447 395 UV filter for TL RC™/RCI™ base
- 10 447 393 ND filter (neutral density filter) for TL RC™/RCI™ base
- 10 450 079 Daylight filter for TL RCI™ base

Controls

- 10 447 443 Leica USB mouse, freely programmable, five-button mouse for connection to TL RCI™ transmitted-light base or PC
- 11 505 180 Leica SmartMove control unit for
- 10 447 444 Leica IsoPro™ motorized mechanical stage
- 10 447 444 Leica PSC control unit for
- 10 450 052 Leica IsoPro™ motorized mechanical stage
- SmartTouch™, external control unit with integr. touchscreen for status control and control of all settings and functions
- 10 450 298 Manual controller for motorized focus
- 33 000 038 USB cables

Ergonomic accessories

- 10 447 431 Leica ErgoRest (handrest for fatigue-free work)



Swiveling arm stands

- 10 447 260 Baseplate, small
- 10 446 436 Baseplate, medium
- 10 447 008 Vertical column 470/35 mm
- 10 447 097 Horizontal arm ESD
- 10 447 098 Horizontal arm standard
- 10 447 006 Flange
- 10 447 007 Stage clamp
- 10 446 437 Baseplate, large
- 10 447 230 Vertical column 800/57 mm
- 10 447 014 Vertical column 560/57 mm
- 10 447 099 Horizontal arm, large
- 10 450 217 Focusing drive with inclinable column
- 10 450 048 Focusing drive coarse/fine for swiveling arm stand



- 10 450 260 Universal plate XL for specimens up to 300 × 300 mm
- 10 450 297 Adapter for universal plate 10 450 260 for all swiveling-arm columns

Leica Z6 APO & Z16 APO

Dimensions of Zoom System, Coaxial Incident Light Housing






Leica Z6 APO









Leica Z16 APO

Coaxial incident light housing










Leica Z6 APO A & Z16 APO A

Dimensions of Zoom System





Dimensions of Objectives, Adapters







Adapters for achromatic objectives

Adapters for DM objectives


Dimensions of Tubes



AS tube





A tube



Y tube


Dimensions of Carriers

Carrier for AS tube




- △ Interface for CMO Focus Drivers
- △ Interface for AS Tube 10447204

Tiltable focusing arm



Focusing arm for probers



“With the user, for the user”

Leica Microsystems

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

• Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

• Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

• Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

• Medical Division

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, “with the user, for the user,” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Active worldwide

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Belgium:	Groot Bijgaarden	Tel. +32 2 790 98 50	Fax +32 2 790 98 68
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Ballerup	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Nanterre Cedex	Tel. +33 811 000 664	Fax +33 1 56 05 23 23
Germany:	Wetzlar	Tel. +49 64 41 29 40 00	Fax +49 64 41 29 41 55
Italy:	Milan	Tel. +39 02 574 861	Fax +39 02 574 03392
Japan:	Tokyo	Tel. +81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives in more than 100 countries

In accordance with the ISO 9001 certificate, Leica Microsystems (Switzerland) Ltd, Industry Division, has at its disposal a management system that meets the requirements of the international standard for quality management. In addition, production meets the requirements of the international standard ISO 14001 for environmental management.