From Eye to Insight





LEICA LASER MICRODISSECTION SYSTEMS

Laser Microdissection (LMD) is a technology for precision sample preparation. In many areas of research it is a basic prerequisite for obtaining well-defined starting material for downstream experiments. Meaningful analyses in the fields of **genomics**, **transcriptomics**, **microarrays**, **next generation sequencing**, **biochips**, **and proteomics** are attained using this high-precision technology.

The Leica LMD systems perform sample preparation for molecular biology analysis directly from the tissue section using a UV laser. The development of innovative new methods and new instruments has made laser microdissection popular in additional fields such as live cell research, climate research, and cover-slip engraving for electron microscopy. Now, more than ever, researchers are using LMD to maximize their research impact.

SUPERIOR SOLUTION

Leica Microsystems offers an extremely precise, highly selective laser microdissection method for a broad range of applications. Users of LMD use our superior systems for various applications:

- Fast, precise isolation of ultra-pure cells and cell populations
- High quality dissectates for genomics, transcriptomics, proteomics, metabolomics, and live cell applications
- Convenient laser manipulation of live cells and other samples
- Mark and track microscopic samples or sample holders

METHOD OF CHOICE

The Leica LMD6 and 7 employ a very gentle, laser-based microscopic sample manipulation, dissection, and collection technique:

- Patented laser beam movement via optics*
- Real-time live cutting directly on the sample with a Pen-screen
- Patented specimen collection via gravity**
- > Dedicated objectives for LMD
- › Adjustable laser settings
- High-end, fully automated upright research microscope
- › Easy-to-use LMD Software
- Patented EP 1276586, US 7035004, JP 3996773
- * Patented DE 10057292, EP 1207392, JP 3641454



The Leica LMD process (from top to bottom):

Step 1: Define the region of interest

Step 2: Laser beam pre-

cisely steered by prisms

along your definition

Step 3: Dissectate is collected by gravity



DISSECTION PERFECTION

Observe and obtain ultrapure and homogenous samples from heterogeneous starting material – contact- and contamination-free simply by gravity. The Leica LMD6/7 systems make this possible: they combine a fully automated upright high-end research microscope and a UV laser.

TWO SYSTEMS, TWO DIFFERENT LASERS

Leica offers two freely configurable laser microdissection systems: the Leica LMD6 (Fig. left) and the Leica LMD7 (Fig. right). Both systems are based on a high-end upright research microscope (Leica DM6 B). In addition to the advanced functions of this fully automated research microscope, Leica LMD systems enable you to isolate and manipulate microscopic samples in a contact- and contamination-free manner. The difference between the two systems is the laser power and flexibility. The laser settings of both systems can be adjusted to perfectly match the needs for your application.

- The LMD6 is ideal for standard laser microdissection applications, e.g. reliable single cell or tumor isolation from soft tissue sections.
- The LMD7 has a high power laser with additional options to adjust the laser pulse frequency and laser head current for dissection of hard tissues like bone, teeth, wood or plant tissue as well as chromosomes.

LASER COMPARISON

LMD6	LMD7
355 nm	349 nm
80 Hz	10-5.000 Hz
< 4 ns	< 4 ns
70 µJ	120 µJ
	355 nm 80 Hz < 4 ns





ADVANTAGES OF LASER MICRODISSECTION

Leica Laser Microdissection -

simple, intuitive, gentle, and smart sample isolation, collection, and manipulation.

> Laser beam movement via optics*:

- Fast, precise, and reliable laser cuts
- Real-time cutting while sample remains fixed
- Convenient documentation by time lapse movies

> Specimen collection via gravity:

- Simple, gentle, contact- and contamination-free
- Allows standard consumables for collection
- No limitation of size or shape of dissectate
- Pool unlimited amounts of dissectates

Adjustable high-powered laser:

- Flexible for a variety of specimens and applications
- Full control of laser power, laser aperture, laser speed, laser frequency (Leica LMD7), and laser focus

> Fully integrated fluorescence:

- Specially designed fluorescence filter cubes
- Live cutting within brightfield and fluorescence

> Upright microscope:

- All benefits of a fully automated high-end research microscope
- Safer and smart dissectate collection by gravity without any additional force
- Free choice between LED or Halogen illumination

Specially designed LMD objectives:

- Ensure the highest possible laser power
- Range of dedicated LMD objectives: 5x, 6.3x, 10x, 20x, 40x, 63x and 150x
- 1.25x for fast slide overviews
- Objectives for additional applications (other than LMD applications, e.g. 100x oil for FISH)

> Easy-to-use software:

- Simple, time-saving, and workflow-based system functionality
- Additional modules for automated pattern recognition (AVC) and automated image capture and documentation (LIF database)

› Leica camera range:

- Choose a camera specific to your application
- Attach up to two cameras (e.g. one for fluorescence, one for brightfield)



* Patented EP 1276586, US 7035004, JP 3996773

** Patented DE 10018255, JP 4236844

CONSUMABLES

Leica Microsystems offers application specific consumables for laser microdissection with different types of membranes on metal frames, glass slides, and ibidi[®] and Petri dishes in different sizes. Whether you need slides that have no autofluorescence, use DIC contrast to view the specimen, or need a suitable surface for growing cell cultures – Leica Microsystems has **the perfect solution for your application**.

FOR DISSECTION WE RECOMMEND

- $\scriptstyle >$ Leica MembraneSlides: glass slides with PEN or PPS membrane
- Leica FrameSlides: steel frames with PEN, PPS, PET, POL or FLUO membrane
- > Leica CoverslipSlides: coverslips with PEN membrane
- > All membranes are highly UV-absorbent
- > DIRECTOR[™] slides
- › Petri dishes with PEN membrane
- › ibidi® slides (e.g. 18-well) with PEN membrane*
- Membrane rings with PEN membrane*
- In addition, Leica LMD Systems support collection from plain glass slides (Draw+Scan ablation)
- * suitable for stack dissection and collection

RANGE OF COLLECTION DEVICES

For motorized and scanning stages:

- > Standard PCR tube caps 0.2 ml
- > Standard PCR tube caps 0.5 ml
- › 8-strip wells
- > Petri dishes with or without PEN membrane

For scanning stage only:

- > 8-strip caps (suitable for 8-strip wells and 96-well plates)
- ibidi[®] slides
- > Multi-well slides (dimensions max. 76 x 26 mm)





STAGES* AND HOLDING DEVICES FOR LMD SYSTEMS



The motorized stage

- > For one slide and Petri dishes
- * Patented DE 10018251, US 6907798, JP 4146642



The scanning stage

- > High flexibility for specimen holders and collection devices
- > High speed and precision

HOLDERS



Holder for 1 slide (25 x 76 mm) 11532732



Holder for big slide (50 x 76 mm) 11505214

Collector with lever

(4 x 0.5 ml PCR tubes)

for tube caps**

11505130



Holder for Petri dish with PEN membrane bottom 11505257



Holder for 3 slides (25 x 76 mm) 11505226



Holder for big slide (50 x 76 mm) and Petri dish with PEN membrane bottom 11505227




Holder for 18-well ibidi® slides and slide stacks 11505255

COLLECTORS




Collector with lever for tube caps** (4 x 0.2 ml PCR tubes) 11505131





Collector for 8-well strip tubes (high volume) 11505258

Collector for tube caps (4 x 0.2 ml PCR tubes) 11505229



Collector for tube caps (4 x 0.5 ml PCR tubes) 11505228




Universal collector for 8-strip tubes, 8-strip caps, slides (Multi-well slides, ibidi®, etc.) and chamber slides 11505276

UNIQUE LASING WITHIN FLUORESCENCE


The ability to perform fluorescence imaging and laser microdissection simultaneously is becoming more important. Lasing within fluorescence using the unique patented axisand filter-system* is one of the strengths of the Leica LMD systems. Whether you are inspecting and dissecting immunolabeled tissue sections or live cells expressing a fluorescent protein – the laser microdissection systems Leica LMD6 and Leica LMD7 make this a standard procedure. In addition, the fully automated fluorescence axis minimizes bleaching effects, accelerates your processes, and offers the exactly reproducible experiment conditions if required.

EXAMPLE


Phloem tissue and stone cells Selection Excision



Dissectate in vessel A



Selection






Dissectate in vessel B Other phloem tissue



LMD FLUORESCENT FILTER SYSTEMS

The Leica's range of specialized LMD fluorescence filter systems is continuously growing. These LMD filter systems are from the ET series, are sputtered, and feature impressively steep edges of the excitation and emission spectrum.



Several newly developed fluorescence filter systems for LMD allow the unique simultaneous cutting under fluorescence:

- › Leica LMD-BGR
- › Leica LMD-GFP band pass
- > Leica LMD-GFP long pass
- Leica LMD-Cy3
 Leica LMD-DAPI
- Leica LMD-GFP/Cy3
 Leica LMD-YFP

> Leica LMD-CFP

> Leica LMD-Alexa594

› Leica LMD-Cy5

* Patented EP 1719998, US 7485878

DEDICATED LASER MICRODISSECTION OBJECTIVES

Leica Microsystems offers a portfolio of dry objectives dedicated for Leica LMD systems. These special LMD objectives feature the highest possible UV-transmission and outstanding imaging performance – the Leica SmartCut series (5x–150x).

Objective	Mag.	NA	WD (mm)	BF, POL	FL	DIC, PH	LMD	CI	SOV
HCX PL FLUOTAR**	1.25x	0.04	3.7	+	+	_	_	_	+++
PLAN**	4x	0.1	26.2	+	+	_	+	++	++
UVI	5x	0.12	11.7	+	+	_	+++	++	++
HI PLAN	6.3x	0.13	12.8	+	÷	_	+++	++	++
HCX PL FLUOTAR	10x	0.3	11.0	+	+	+	+++	+	+
UVI	10x	0.25	2.9	+	+	+	+++	-	+
HCX PL FLUOTAR	20x	0.4	6.9	+	+	+	+++	+	+
HCX PL FLUOTAR	40x	0.6	3.3–1.9	+	+	+	+++	_	_
HCX PL FLUOTAR	63x	0.7	2.6–1.8	+	+	+	+++	_	_
HCX PL FLUOTAR	150x	0.9	0.25	+	+	+	+++	_	_

Mag. – magnification; NA – numerical aperture; WD – working distance; BF – brightfield; POL – polarized light; FL – fluorescence; DIC – differential interference contrast; PH – phase contrast; LMD – laser microdissection; CI – cap inspection; SOV – specimen overview; + = suitable; – = not suitable; ++ = dedicated for; +++ = most suitable

** Additionally recommended objectives

OTHER OBJECTIVES

In addition to the LMD objectives, any other Leica objectives can be used with the Leica DM6 B upright microscope for specific applications (e.g. for FISH). In the table above, the 1.25x objective is used for fast specimen overviews and the 4x objective is dedicated for cap inspection.



CAMERA PORTFOLIO

Transmitted light, fluorescence or both: Leica LMD systems and software support a range of digital cameras for different requirements.

Choose the new Leica LMD CC7000 digital color camera dedicated for LMD applications. The highly compact LMD color camera Leica LMD CC7000 (# 11501478) is a GigE camera with 1/3" interline progressive scan CCD sensor and 1.2 Megapixel resolution. Experience the speed of the new Leica LMD CC7000, unique for Leica LMD systems.

Leica LMD CC7000 digital color camera



Choose from the portfolio of Leica microscope cameras



Left: Leica LMD6 with one camera: The Leica LMD CC7000 for ultra fast digital live images

Right: Leica LMD7 with two cameras: The Leica DFC7000 T for ultra fast digital live images and the Leica DFC365 FX for demanding fluorescence.

Dual camera support allows you to combine any two cameras, thus making your system a multifunctional tool for any kind of application!

For Leica LMD6 the set-up with 2 cameras is different.






LIVE CELL CUTTING ACCESSORIES


The Leica LMD systems are ideal for any kind of live cell application. With our "sandwich-technology" using membrane rings or our ibidi[®] slides' sterile microchambers^{*}, live cells are easily dissected in a sterile environment. The LMD system can also be equipped with a climate chamber and full climate control. The time lapse movie function in combination with the laser beam movement by prisms makes it the perfect tool for your laser manipulation experiments. Collect cells or cell clusters via gravity, without any additional force or handling steps, directly into culture media for immediate recultivation, or, alternatively, into collectors for downstream analysis.

APPLICATION EXAMPLE: TRANSGENIC FLUORESCENT CELLS


Human forensic fibroblasts infected with human cytomegalo virus, HCMV-GFP fusion protein (Courtesy of Margarete Digel and Dr. Christian Sinzger, Institute of Medical Virology, UKT University of Tübingen, Germany)




Before microdissection




After microdissection



Infection of microdissected fibroblasts



4 days after recultivation



Leica LMD climate chamber with environmental control

Classical selection of infected cells followed by dilution series takes 2 months.

Patented** universal holder with different slides suitable for live cell applications




* Patented US 7807108

** Patented DE 102009029078

EASY-TO-USE SOFTWARE

Leica LMD Software offers the functions necessary for perfect Laser Microdissection or Laser Manipulation. With full control of all laser settings you can adjust the laser to your specimen, independent of shape or size.

Improve your experiments with the different drawing tools and cutting modes. Enjoy the intuitive user interface and the speed of the system and fully automated fluorescence, DIC, PH and POL.




FEATURES ARE ALREADY INCLUDED IN THE LMD CORE SOFTWARE, FOR EXAMPLE:

- Full microscope control including illumination method, camera control, specimen and collector holders
- Full control of all laser settings for adjustment to any type of sample
- Specimen overview images for optimal orientation and navigation
- Time lapse movie function for application recordings, the sample stays fixed even during laser applications
- Different cutting modes and drawing tools to either dissect and pool or separate specimens of different sizes and shapes
- Annotation and length measure tool for documentation
- Image capture function for documentation
- Save and restore application settings to ensure quick and easy experiment starts

SOFTWARE MODULES


RETAIN ALL INFORMATION: LEICA LMD DATABASE

- Automated storage of images prior to cut, after cut, and after cap inspection
- > Full access to all experiment parameters
- Shared database LIF-file format for image access and processing with Leica LAS AF (Leica Application Suite, Advanced Fluorescence)
- Image export to current image file formats (JPG, TIF, ...)




SPEED UP YOUR LMD APPLICATION: AVC (AUTOMATED CELL RECOGNITION*) FOR PATTERN RECOGNITION

- Automated marking of cutting lines to avoid time-consuming manual selection
- Fully integrated solution for automated pattern recognition
- > Easy-to-use interface
- Works well with any kind of specimen (transmitted light and fluorescence)
- Option to save and restore different settings
- Smart autofocus for best performance in different field of views
- * Patented EP 1676116



SUMMARY Maximize your impact with a Leica LMD system



MAIN APPLICATION FIELDS FOR LEICA LASER MICRODISSECTION SYSTEMS


- › Cancer research
- › Pathology
- > Molecular biology research
- > Neuroscience research
- > Developmental research
- > Plant research
- › Forensic research
- > Physiology research
- Clinical research
- > Pharmaceutical research
- > ... almost anywhere

MAIN APPLICATIONS FOR LEICA LASER MICRODISSECTION SYSTEMS

- Extraction of homogeneous samples from heterogeneous starting material for genomic, transcriptomic, proteomic, and metabolite analysis
- Live cell cloning, manipulation, and re-culturing
- Ablation and damage of live cells, tissues, and embryos monitored by time-lapse movies
- > Thrombosis inducement
- TEM sample selection before resin embedding
- Engraving coverslips for CLEM application
- NanoSIMS support
- > ... your own unique application

BENEFIT FROM MORE THAN 10 YEARS OF EXPERIENCE

- Third generation of proven Leica LMD systems
- Workflow-based, easy-to-use, powerful software
- Large, fast growing library of scientific publications
- Knowledgeable, experienced Leica support personnel





with the market leaders.

LIFE SCIENCE DIVISION

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

INDUSTRY DIVISION

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

MEDICAL DIVISION

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

Leica Microsystems operates globally in three divisions, where we rank Leica Microsystems - an international company with a strong network of worldwide customer services:

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Diegem	+32	2 790 98 50	2 790 98 68
Brazil · São Paulo	+55	11 2764-2411	11 2764-2400
Canada · Concord/Ontario	+1	800 248 0123	847 405 0164
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
India · Mumbai	+91	226 1880 200	226 1880 333
ltaly · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 6758 5630	3 5155 4333
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People's Rep. of China · Hong Kong	+852	2564 6699	2564 4163
· Shanghai	+86	21 6039 6000	21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6550 5999	6564 5955
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Bromma	+46	8 6 2 5 4 5 4 5	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
Turkey · Istanbul	+90	216 504 0100	216 504 0110
United Kingdom · Milton Keynes	+44	800 298 2344	1908 577640
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 405 0164



http://www.leica-microsystems.com/products/ light-microscopes/life-science-research/ laser-microdissection/

The Leica LMD6/7 microscopes comply with the Council Directive 98/79/EEC concerning in vitro diagnostics.

They also conform to the Council Directives 2006/95/EG concerning electrical apparatus and 2004/108/EG concerning electro-magnetic compatibility for use in an industrial environment.