



CERTIFICATE OF ANALYSIS

Product **R.T.U HORSERADISH PEROXIDASE AVIDIN D**

Catalog No. A-2704

Amount 100 ml

Lot No. ZC0223

Concentration 1 µg/ml active conjugate*

Molar ratio of active conjugate

$$\frac{\text{moles of HRP}}{\text{moles of Avidin}} = 1.5$$

Storage Conditions Refrigerate; DO NOT FREEZE

Remarks:

* This product is provided at a ready-to-use concentration in a solution which preserves the activity on the conjugate when stored as recommended.

For most applications the recommended incubation time is 30 minutes.

Do not use sodium azide with this product since it is an inhibitor of horseradish peroxidase.

FOR LABORATORY USE.