



---

## CERTIFICATE OF ANALYSIS

Product AGAROSE DATURA STRAMONIUM LECTIN

Catalog AL-1183

Amount 2 ml of settled gel (1:1, v/v slurry)

Lot No. W0203

Suspension Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.1 mM Ca<sup>++</sup>,  
0.08% sodium azide

Protein Concentration 3.0 mg/ml of settled gel

Binding Capacity 1.6 mg asialo-fetuin/ml of gel

Storage Conditions Refrigerate - DO NOT FREEZE

Other Remarks

Chitin hydrolysate (Cat. No. SP-0090) is recommended for elution.