



CERTIFICATE OF ANALYSIS

Product **AGAROSE GRIFFONIA (BANDIERAEA) SIMPLICIFOLIA LECTIN II**

Catalog AL-1213

Amount 2 ml of settled gel (1:1, v/v slurry)

Lot No. V0127

Suspension Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 20 mM GlcNAc,
0.08% sodium azide

Protein Concentration 3.0 mg/ml of settled gel

Binding Capacity Greater than 4.0 mg of N-acetylglucosaminyl glycoprotein/ml gel

Storage Conditions Refrigerate - DO NOT FREEZE

Other Remarks

Wash gel thoroughly with buffer before use to remove sugar added to stabilize the lectin.

Recommended conditions to elute bound glycoconjugates: 0.2M N-acetylglucosamine in buffered saline or Chitin Hydrolysate (Cat. No. SP-0090).