



---

## CERTIFICATE OF ANALYSIS

Product AGAROSE AGARICUS BISPORUS LECTIN  
(White Button Mushroom) (ABL)

Catalog No. AL-1423

Amount 2 ml of settled gel (1:1, v/v slurry)

Lot No. ZB0810

Suspension Solution 10 mM HEPES, pH 7.5, 150 mM NaCl, 0.08% sodium azide

Protein Concentration 2 mg/ml of settled gel

Binding Capacity > 1.5 mg fetuin/ml gel

Storage Refrigerate - DO NOT FREEZE

Purity Made from affinity purified lectin which is homogeneous by SDS-PAGE.

### Other Remarks

Sugar Specificity: Gal( $\beta$ -1,3)GalNAc (T disaccharide)

Elution conditions: Fetal calf fetuin bound to Agarose *Agaricus bisporus* Lectin can be eluted with 0.1 M ammonium hydroxide. However, because glycoproteins vary in their structure, some glycoproteins may not elute under these conditions.