



CERTIFICATE OF ANALYSIS

Product ANTIBODY TO SOYBEAN AGGLUTININ
 AFFINITY PURIFIED

Catalog AS-2014

Amount 1 mg

Concentration 1mg/ml

Lot No. ZC0825

Produced in goat

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has
 been added for stabilization.

Storage Refrigerate; for long term storage, aliquots may be stored frozen

Active antibody by solid-phase binding assay 99 %

Remarks

The recommended concentration range for use is 1 to 10 µg/ml.