

CERTIFICATE OF ANALYSIS

Product ANTIBODY TO WHEAT GERM AGGLUTININ

AFFINITY PURIFIED

Catalog No. AS-2024

Amount 1 mg

Concentration 2mg/ml

Lot No. <u>ZC0621</u>

Produced in goat

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has

been added for stabilization.

Storage Refrigerate: for long term storage, aliquot and store frozen

Active antibody by solid-phase binding assay 99 %

Remarks

The recommended concentration range for use is 1 to $10 \mu g/ml$.