



---

## CERTIFICATE OF ANALYSIS

Product      **ANTIBODY TO PEANUT AGGLUTININ**  
                 **\*\*AFFINITY PURIFIED\*\***

Catalog      AS-2074

Amount      1 mg

Concentration   2 mg/ml

Lot No.        ZC1207

Produced in    goat

Solution       10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage        Refrigerate; for long term storage, aliquot and store frozen

Active Antibody by solid-phase binding assay   99%

Remarks:

The recommended concentration range for use is 1 to 10 µg/ml.