



CERTIFICATE OF ANALYSIS

Product **ANTIBODY TO *GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA* LECTIN I**
****AFFINITY PURIFIED****

Catalog AS-2104

Amount 1 mg active conjugate

Concentration 2 mg/ml

Lot No. Z0724

Produced in goat

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08%
sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long term storage, aliquot and store frozen

Active antibody by solid-phase binding assay 91%

Remarks:

The recommended concentration range for use is 1 to 10 $\mu\text{g/ml}$.

When reconstituting, do not vortex.