



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED WHEAT GERM AGGLUTININ**

Catalog B-1025

Amount 5 mg active conjugate

Lot No. Z0611

Concentration 5 mg/ml

Active Conjugate 98%

Lectin Biotinylated 99%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen.
If a precipitate forms upon long-term storage at 4 °C, warm to
37 °C before use.

Solution: 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca⁺⁺.
For most applications we recommend a freshly prepared working solution of 5-20 µg/n
in the above buffer.