



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *LENS CULINARIS* AGGLUTININ**

Catalog B-1045

Amount 5 mg

Concentration 5 mg/ml

Lot No. ZC1221

Active Conjugate 96%

Lectin Biotinylated 98%

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.01 mM Mn⁺⁺, 0.1 mM Ca⁺⁺, 0.08% sodium azide.

Storage Conditions Refrigerate. Lectin solutions may be aliquoted and stored frozen or refrigerated.

For most applications we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.