



---

## CERTIFICATE OF ANALYSIS

Product     **BIOTINYLATED *PISUM SATIVUM* AGGLUTININ**

Catalog     B-1055

Amount     5 mg

Concentration     5 mg/ml

Lot No.     Y0104

Active Conjugate     99%

Lectin Biotinylated     99%

Solution     10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.1 mM Ca<sup>++</sup>, 0.01 mM Mn<sup>++</sup>,  
0.08% sodium azide

Storage Conditions     Refrigerate; for long term storage, aliquot and store frozen.

For most applications we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.