



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED PEANUT AGGLUTININ**

Catalog No. B-1075

Amount 5 mg

Concentration 5 mg/ml

Lot No. X1221

Active Conjugate 98%

Lectin Biotinylated 99%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺,
0.01 mM Mn⁺⁺

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen.

Remarks:

DO NOT VORTEX.

For most applications we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer.

Use of buffers containing 0.1 mM Ca⁺⁺ and 0.01 mM Mn⁺⁺ is recommended.