



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *RICINUS COMMUNIS* AGGLUTININ I (RCA₁₂₀)**

Catalog **B-1085**

Amount **5 mg**

Lot No. **Y0214**

Active Conjugate **96%**

Lectin Biotinylated **99%**

Protein Concentration **5.0 mg/ml**

Solution **10 mM phosphate, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.**

Storage **Refrigerate; DO NOT FREEZE**

Remarks:

EXERCISE CAUTION IN ITS USE! It has been found to have moderate toxicity in some animal models, but its toxicity to humans has not been fully determined, and a related protein, RCA 60, is highly toxic to humans. Never pipet solutions of this lectin by mouth.
DO NOT LYOPHILIZE!

The recommended concentration range for use is 5-20 μ g/ml.