



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED PHASEOLUS VULGARIS LEUCOAGGLUTININ (PHA-L)**

Catalog B-1115

Amount 2 mg active conjugate

Lot No. X0804

Lectin Biotinylated 99%

Concentration 2.0 mg/ml

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺

Storage Refrigerate

Other Remarks:

For most applications we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer.