



## CERTIFICATE OF ANALYSIS

Product     **BIOTINYLATED *SOLANUM TUBEROSUM* (POTATO) LECTIN**

Catalog     B-1165

Amount     2 mg

Lot No.     X0502

Active Conjugate   91%

Lectin Biotinylated   99%

Solution     10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca<sup>++</sup>.

Storage Conditions   Refrigerate; for long term storage, aliquots may be stored frozen

Other Remarks:

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.