



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *DATURA STRAMONIUM* LECTIN**

Catalog B-1185

Amount 2 mg active conjugate

Concentration 2mg/ml

Lot No. Z1119

Active Conjugate 97%

Lectin Biotinylated 96%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺.

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Other Remarks

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.

*Aggregation may occur with time if the biotinylated lectin is concentrated to greater than 2 mg/ml.