



---

## CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA* LECTIN I  
ISOLECTIN B4**

Catalog B-1205

Amount 0.5 mg active conjugate

Concentration 0.5mg/ml

Lot No. X0819

Active Conjugate 95%

Lectin Biotinylated 99%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Other Remarks:

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.

This lot will agglutinate blood group B glutaraldehyde-fixed human red blood cells at a concentration of 3 µg/ml. No agglutination of blood group A, glutaraldehyde-fixed, human red blood cells occurred at concentrations up to 100 µg/ml.