



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA* LECTIN II**

Catalog B-1215

Amount 2 mg

Lot No. X0301

Active Conjugate 99%

Lectin Biotinylated 99%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Other Remarks:

Reconstitute biotinylated lectin by addition 1ml of water. The resulting solution will have the following composition: 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺. For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.