



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *VICIA VILLOSA* LECTIN**

Catalog B-1235

Amount 2 mg active conjugate

Concentration 2 mg/ml

Lot No. ZA0708

Active Conjugate 99%

Lectin Biotinylated 99%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide,
0.1 mM Ca⁺⁺, 0.01 mM Mn⁺⁺.

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.