



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED GALANTHUS NIVALIS (SNOWDROP) LECTIN**

Catalog B-1245

Amount 2 mg

Concentration 2 mg/ml

Lot No. ZA0715

Active Conjugate 98%

Lectin Biotinylated 97%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺.

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage, warm to 37 °C and centrifuge before use.