



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED MAACKIA AMURENSIS LECTIN II (MAL II)**

Catalog B-1265

Amount 1 mg active conjugate

Lot No. Z0611

Concentration 1 mg/ml

Lectin Biotinylated 97%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide. For most applications, we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer.

Storage Conditions Store 4 °C