



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *SAMBUCUS NIGRA* (ELDERBERRY) BARK LECTIN**

Catalog B-1305

Amount 2 mg active conjugate

Concentration 2 mg/ml

Lot No. Z1002

Active Conjugate 93%

Lectin Biotinylated 99%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca⁺⁺.

Other Remarks:

For most applications, we recommend a freshly prepared working solution of 10 μ g/ml in the above buffer.

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage, warm to 37 °C and centrifuge before use.