



## CERTIFICATE OF ANALYSIS

Product     **BIOTINYLATED LOTUS TETRAGONOLOBUS LECTIN**  
                  from *Tetragonolobus purpureus* seeds

Catalog     B-1325

Amount     2 mg active conjugate

Lot No.     ZC2428

Active Conjugate   98%

Lectin Biotinylated   99%

Storage Conditions   Refrigerate; for long term storage, aliquot and store frozen

Reconstitute by adding 1ml water. The resulting solution will have the following composition: 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca<sup>++</sup>. For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.