



---

## CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED MACLURA POMIFERA LECTIN**

Catalog B-1345

Amount 2 mg active conjugate

Concentration 2mg/ml

Lot No. ZB0211

Active Conjugate 99%

Lectin Biotinylated 99%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide.

Storage Conditions Refrigerate; for long term storage, aliquot and store frozen

Other Remarks:

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.

Inhibiting Sugar: 500 mM galactose (Cat. No. S-9003).

For immunofluorescence applications, we recommend that sections be mounted in VECTASHIELD® Hard+Set™ Mounting Medium (Cat. No. H-1400) or VECTASHIELD® Hard+Set™ Mounting Medium with DAPI (Cat. No. H-1500).