



---

## CERTIFICATE OF ANALYSIS

Product     **BIOTINYLATED *WISTERIA FLORIBUNDA* LECTIN**

Catalog     B-1355

Amount     2 mg

Concentration   2 mg/ml

Lot No.     Y0409

Active Conjugate   98%

Lectin Biotinylated   99%

Solution     10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca++, 0.1 mM Mn++.

Storage Conditions   Refrigerate; for long term storage, aliquot and store frozen

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.