



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *NARCISSUS PSEUDONARCISSUS* (DAFFODIL) LECTIN**

Catalog No. B-1375

Amount 2 mg

Concentration 2 mg/ml

Lot No. ZC0315

Active Conjugate 95%

Lectin Biotinylated 99%

Storage Conditions Refrigerate

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide

Other Remarks

For most applications, we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.