

CERTIFICATE OF ANALYSIS

AMCA ANTI-HUMAN LAMBDA CHAIN **Product**

AFFINITY PURIFIED

Catalog No. CI-3070

0.5 mg active conjugate Amount

ZB0702 Lot No.

Produced in goat

1.0 mg/ml active conjugate Concentration

10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide. Solution

None Free-Fluorochrome

<u>>8</u> F/P (molar)

human lambda IgG

Maximum Excitation @ 345-355 nm Maximum Emission @ 448-454 nm

Refrigerate in the dark. If a precipitate forms, warm to 37 °C before use. Storage

Antibody Activity (by solid phase IgG binding assay) 95%

Solid-phase enzyme immunoassay of cross-reactivity:

Protein Percent Cross-reactivity 50%

human kappa IgG <1%

<1% free kappa chain

If this fluorochrome-labeled antibody is to be used in tissues which may contain cross-reacting endogenous immunoglobulins, dilution of this antibody may be made in buffers containing 2% normal serum from the same species as the tissue.