

Consumables

Chrome rod 2 x 25mm (5 pcs.)

16771790

- For high resolution chrome coatings Dimensions: Ø 2 x 25 mm Pack of 5

Consumables for Carbon Evaporation

16771791

- Consisting of:
 - 50 tungsten-cathodes
 - 30 carbon rods

Consumables Ta/Tu

16771792

- Consumables Ta/Tu

Consumables platinum/carbon e-beam evap

16771793

- Consists of:
50 W-cathodes, 60 Pt inserts, 30 carbon rods (2 sides each usable)

Tungston rod (10 pieces)

16771795

- Tungston rod (10 pieces)

Spectrographically Pure Carbon Rod

16771796

- To produce Carbon/Platinum mix coatings. Both sides drilled for platinum inserts. Dimensions: Ø 2 x 20 mm Pack of 10

Spectrographically Pure Carbon Rod

16771797

- To produce carbon backing films. Dimensions: Ø 3 x 50 mm Pack of 6

Platinum Insert (10p)

16771798

- For Carbon/Platinum mix coatings. Dimensions: Ø 1.5 x 2 mm Pack of 10

Tungsten cathode (5 pieces)

16771799

- suitable for evaporation sources, Pack of 5