CRUMAFL2

Cruma breaks with the tradition of removable laminar flow cabinets offering the new Cruma FL-2, an innovative cabin whole body.

Furthermore, the new cabin Cruma FL-2 combines a striking and unique exterior design with the best technology and internal development. However, to meet the operator in his everyday work and protect it in their work caring for the environment.



USES

- √ Laboratory food industries in general
- √ Transfer of drugs in the pharmacy department
- √ Hematology and clinical analysis
- √ Filled with antibiotics and injectable drugs (except cytostatic) microscopic analysis
- √ Tissue culture bags reconstitution of parenteral methods of assisted fertilization techniques
- √ Micropropagation techniques. Plant cell cultures.
- √ Water quality control
- √ Manufacturing of electronic devices...

NEW FEATURES

More information on the new LCD display



- √ New size 127x34mm display
- √ Type of filter installed, working hours, expiration date and next revision date
- √ Countdown timer
- √ Clock and calendar

New features and components

- √ Initial cycle flow adequacy and final purge cycle
- √ Fault LED
- √ Control of air flow through Microprocessor
- √ HEPA filters with electronic chip
- √ LED illumination

New alarms and scheduled warnings

- √ Next validation warning
- √ Few hours of filter life warning
- √ Countdown timer warning
- √ Expired filter alarm by hours
- √ Expired filter alarm by date
- √ Equipment without filter alarm



MAIN FEATURES

- Manufacture: outer surfaces made of cold rolled steel gray / ocher coating, stoveenamelled, resistant to the chemical compounds used for sterilization
- √ Worktop: the surface of the working chamber is in secured safety glass
- √ Absolute HEPA H-14 filter with flow equalizer, efficiency greater than 99.995% on 0.3 micron particles. Filter is in accordance with EN-779:1996 and EN-1822:1998
- $\sqrt{\mbox{ Fan:}}$ centrifugal double suction, suitable to provide laminar flow in the workspace
- √ Control panel with soft-touch keyboard and switch for LED and optional UV lamps; the air speed is kept at a constant rate thanks to the HEPA filter which is electronically self- regulated, electronic control air speed and control key to slow it down
- $\sqrt{\,$ Lighting: LED lamp placed in the central part of the cabinet, for a better light distribution over the working surface

TECHNICAL INFORMATION		
Number of filtering columns		1
Number of IP44 fans		1
Average face velocity		0,40 m/s
Total electrical power consumption		226 W
Voltage/Frequency		110-220 V / 50-60 Hz
Illumination level LED		900 Lux
UV germicidal lamp		15W
Noise level		50 dB
Packaging: 100% recyclable Wood box with international phytosanitary certificate	Volume	0,95 m ³
	Weight	142 Kg

It is not a typographical error, 7 year warranty

Because we are convinced of the quality of our products.



*More information www.cruma.es

Well done, well shipped. Our responsible packaging

Wood box 100% recyclable with international phytosanitary certificate.



Do you need help or technical assistance?

Contact your distributor of call us if you have any questions or need technical support, spare parts, maintenance service...

+34 93 370 61 62





CRUMAFLE

SIZES (mm) *CRUMAFE2 463 601

SIZES (mm)					
	Exterior			Interior	
Width 1200	Depth 600	Height 1137	Width 1110	Depth 556	Height 610

STRUCTURE	
Metallic parts	1.2mm galvanized steel, coated with antiacid polymer powder resin thermo-hardened at 200 $^{\circ}$ C
Front Panel	Transparent polymethylmethacrylate 6 mm thick (light transmission of 93%)







STANDARD EQUIPMENT		
Electronic circuit with large format LCD display	Security levels: level 1 for users and level 2 for maintenance users	
Filters controlled by microprocessor	The filters are provided with a microchip with mini USB connection identifying the type of filter installed, the expiry date and the serial number	
Illumination	96 LED high light intensity and low consumption 16 Watts / 900 Lux	
Digital timer with audible warning	Hourmeter for the control and warning of the work in the cabinet	
Electrical plug	Placed in the internal roof of the cabinet and conmuted with the electronic circuit	
Clock and Calendar	Display of date and time	
Gastap	Placed in the internal of the cabinet	
Work surface	Liquid retention tray (2-10 liters) with working surface in white tempered glass	
G4 prefilter	Pre-filtering class G4 biofibre synthetic blanket	
Warranty	7 years	

OPTIONAL EQUIPMENT
Tubular stand
Stand with wheels Movilair
Double gas tap
Vacuum tap
Double vacuum tap
Double electric socket

ACCORDING TO STANDARDS		
Filters	EN-779: 1996 (HEPA & ULPA Filters) EN-1822:1998 (HEPA & ULPA Filters)	
Quality	UNE EN ISO 9001:2015	

We recognise our responsability and dependence towards a healthy environment and, therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

our 72 engagement



















DISTRIBUTED BY: