



CERTIFICATE OF ANALYSIS

Product **Cy[®]3 GOAT ANTI-RABBIT IgG (H+L)**
****AFFINITY PURIFIED****

Catalog No. CY-1300

Amount 1 mg active conjugate

Lot No. ZC1208

Concentration 1.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 550 nm
Maximum Emission @ 570 nm


F/P (molar) 4.7

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solid-phase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

human, mouse, horse, cow, pig, cat, dog, chicken, monkey, guinea pig, hamster, rat, sheep.



Cy[®] is a registered trademark of GE Healthcare