



## CERTIFICATE OF ANALYSIS

Product **Cy<sup>®</sup>5 HORSE ANTI-MOUSE IgG (H+L)**  
**\*\*AFFINITY PURIFIED\*\***

Catalog No. CY-2500

Amount 1 mg active conjugate

Lot No. ZC1208

Concentration 1.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 647 nm  
Maximum Emission @ 663 nm


F/P (molar) 4.8

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solid-phase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

**human, rabbit, goat, cow, pig, cat, dog, chicken, monkey, guinea pig.**



*Cy<sup>®</sup> is a registered trademark of GE Healthcare*