



DAW50

DAI recognized leadership in **full-plate Microplate washers** extends to the DAW50 Auto Strip Washer. The DAW50 compact footprint contains a powerhouse of washing capabilities unsurpassed in its class. The washer's excellent dispense accuracy and evacuation efficiency supports 96- and 384-well strip or plate washing. Comprehensive and easy-to-use onboard software provides the utmost flexibility for a variety of applications. CE certified.

Microplate Types:

All flat-bottom 96-well plates, 384-well flat-bottom plates 1 x 8 flat-bottom strips (16-channel), All 96-well plates, 1 x 8 strips (8-channel)

All 96-well plates, 1 x 12 strips (12-channel)

Mainfold Type:

16-channel patented Dual-Action, with separate dispense/aspiration (16-well), 8-channel aspiration/dispense (8-channel), 12-channel aspiration/dispense (12-channel)

Washing Speed: 3 asp./disp. Cycles:

< 250 second (100 μ l/well, 384 wells)

> 105 second (300 μ l/well, 96 wells)

Volume Range: 25 to 3000 μ l, Soak Time: 1 to 600 seconds

Power: 100 to 240 VAC, 50 to 60 H