



## CERTIFICATE OF ANALYSIS

Product DyLight<sup>®</sup> 488 ANTI-RABBIT IgG (H+L)  
\*\*AFFINITY PURIFIED\*\*

Catalog No. DI-1088

Amount 1.5 mg active conjugate

Lot No. ZC1005

Produced in horse

Concentration 1.5 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 493 nm

Maximum Emission @ 518 nm


F/P (molar) 6.7

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solid-phase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

**human, mouse, goat, pig, cat, dog, chicken, monkey, guinea pig, hamster, rat, sheep.**



*DyLight<sup>®</sup> is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.*