

CERTIFICATE OF ANALYSIS

Product DyLight[®]549 ANTI-RABBIT IgG (H+L) <u>**AFFINITY PURIFIED**</u>

Catalog No. DI-1549

Amount <u>1.5 mg active conjugate</u>

Lot No. <u>ZB0506</u>

Produced in goat

Concentration <u>1.5 mg/ml active conjugate</u>

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 556 nm

Maximum Emission @ 571 nm


F/P (molar) <u>7.5</u>

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solidphase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

human, mouse, horse, cow, pig, cat, dog, chicken, monkey, guinea pig, hamster, rat, sheep.



DyLight[™] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Vector Laboratories, Inc., 30 Ingold Road, Burlingame, CA 94010 U.S.A. Tel: (650)697-3600 • Fax: (650)697-0339 • Email: vector@vectorlabs.com • Web: www.vectorlabs.com