



CERTIFICATE OF ANALYSIS

Product **DyLight[®] 594 ANTI-RABBIT IgG (H+L)**
****AFFINITY PURIFIED****

Catalog No. DI-1594

Amount 1.5 mg active conjugate

Lot No. ZA0618

Produced in goat

Concentration 1.5 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 592 nm

Maximum Emission @ 617 nm

F/P (molar) 4.9


Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solid-phase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

**human, mouse, horse, cow, pig, cat, dog, chicken, monkey,
guinea pig hamster rat sheen**

DyLight[®] 594 Conjugate Spectral Properties



DyLight[®] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.