

CERTIFICATE OF ANALYSIS

DyLight®649 ANTI-RABBIT IgG (H+L) Product

AFFINITY PURIFIED

DI-1649 Catalog No.

1.5 mg active conjugate Amount

Z0820 Lot No.

Produced in goat

1.5 mg/ml active conjugate Concentration

10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide. Solution

Maximum Excitation @ 655 nm Maximum Emission @ 670 nm


<u>4.7</u> F/P (molar)

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solidphase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

human, mouse, horse, cow, pig, cat, dog, chicken, monkey, guinea pig, hamster, rat, sheep.



DyLight[®] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.