

CERTIFICATE OF ANALYSIS

Product <u>DyLight®549 ANTI-MOUSE IgG (H+L)</u>

AFFINITY PURIFIED

Catalog No. DI-2549

Amount <u>1.5 mg active conjugate</u>

Lot No. <u>Y1030</u>

Produced in horse

Concentration 1.5 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 556 nm Maximum Emission @ 571 nm


F/P (molar) 6.3

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-20 µg/ml

The affinity purified antibody has been extensively adsorbed on solid phases and solidphase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

human, rabbit, goat, cow, pig, cat, dog, chicken, monkey, guinea pig.



DyLight[®] *is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.*