

CERTIFICATE OF ANALYSIS

Product <u>DyLight®649 ANTI-MOUSE IgG (H+L)</u>

AFFINITY PURIFIED

Catalog No. <u>DI-2649</u>

Amount <u>1.5 mg active conjugate</u>

Lot No. <u>ZA0424</u>

Produced in horse

Concentration 1.5 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 655 nm Maximum Emission @ 670 nm

F/P (molar) 4.6

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use $5-20 \mu g/ml$

The affinity purified antibody has been extensively adsorbed on solid phases and solidphase enzyme immunoassay shows minimal cross-reactivity against the following immunoglobulins:

human, rabbit, goat, cow, pig, cat, dog, chicken, monkey, guinea pig.

