



CERTIFICATE OF ANALYSIS

Product **DyLight® 594 ULEX EUROPAEUS AGGLUTININ I**
****AFFINITY PURIFIED****

Catalog No. **DL-106Z**

Amount **1 mg active conjugate**

Lot No. **ZC1205**

Concentration **1 mg/ml**

Solution **10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,**
0.1 mM Ca++.

Maximum Excitation @ 592 nm
Maximum Emission @ 617 nm

F/P (molar) **1.7**


Active Conjugate **98%**

Storage Conditions **Refrigerate in the dark.**

Other Remarks:

If a precipitate forms upon long-term storage, warm to 37 °C.

The recommended concentration range for use is 5-20 µg/ml.



DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.