

## **CERTIFICATE OF ANALYSIS**

## Product FLUORESCEIN ANTI-HUMAN IgG (GAMMA CHAIN SPECIFIC) \*\*AFFINITY PURIFIED\*\*

Catalog No. FI-3080

Amount <u>0.5 mg active conjugate</u>

Lot No. <u>W0310</u>

Produced in goat

Concentration <u>1.0 mg/ml active conjugate</u>

Solution <u>10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.</u>

Free-Fluorochrome None

Maximum Excitation @ 490-500 nm Maximum Emission @ 510-520 nm

F/P (molar) 5.4

Storage Conditions <u>Refrigerate in the dark.</u>

Antibody Activity (by solid phase IgG binding assay) 95%

Solid-phase enzyme immunoassay of cross-reactivity:

<u>Protein</u>human IgA, human kappa chain, human IgM, human IgE, human lambda chain

Percent Cross-reactivity

<1%

If this fluorescein-labeled antibody is to be used in tissues which may contain cross-reacting endogenous immunoglobulins, dilution of this antibody may be made in buffers containing 2% normal serum from the same species as the tissue.

The recommended concentration range for use is  $5-20 \ \mu g/ml$ .