



---

## CERTIFICATE OF ANALYSIS

Product FLUORESCHEIN CONCAVALIN A

Catalog No. FL-1001

Amount 25 mg

Lot No. W0909

Concentration 5 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.1 mM Ca<sup>++</sup>, 0.08% sodium azide, 0.01 mM Mn<sup>++</sup>.

Free Fluorescein None

F/P (molar) 5.0

Active Conjugate 99%

Storage Conditions Refrigerate in the dark.

### Other Remarks:

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage at 4 °C, warm to 37 °C and centrifuge before use.

The recommended concentration range for use is 5-20 µg/ml.