



CERTIFICATE OF ANALYSIS

Product FLUORESCIN *GALANTHUS NIVALIS* (SNOWDROP) LECTIN

Catalog No. FL-1241

Amount 2 mg active conjugate

Lot No. Z0412

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 6.0

Active Conjugate 99%

Storage Conditions Refrigerate in the dark

Other Remarks:

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage, warm to 37 °C and centrifuge before use.

The recommended concentration range for use is 5-20 µg/ml.