



CERTIFICATE OF ANALYSIS

Product **FLUORESCCEIN *AGARICUS BISPORUS* LECTIN
(White Button Mushroom) (ABL)**

Catalog No. FL-1421

Amount 1 mg

Concentration 1 mg/ml

Lot No. ZB0811

Solution 10 mM HEPES, pH 7.5, 150 mM NaCl, 0.08% sodium azide

Storage Refrigerate

Purity Made from affinity purified lectin which is homogeneous by SDS-PAGE.

Other Remarks:

Sugar Specificity: Gal(β -1,3)GalNAc (T disaccharide)

For immunofluorescence applications, we recommend that tissue sections be mounted in VECTASHIELD[®] Hard•Set[™] Mounting Medium (H-1400) or VECTASHIELD[®] Hard•Set[™] Mounting Medium with DAPI (H-1500).