



---

## CERTIFICATE OF ANALYSIS

Product ANTI-HUMAN IgG (gamma chain specific)  
\*\*AFFINITY PURIFIED\*\*

Catalog AI-3080

Amount 0.5 mg IgG

Concentration 1 mg/ml

Lot No. ZB0515

Produced in goat

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long-term storage, aliquot and store frozen.

Antibody Activity (by solid phase IgG binding assay) 91%

Solid-phase enzyme immunoassay of cross-reactivity:

<u>Protein</u>	<u>Percent Cross-reactivity</u>
human IgE, human IgM	<1%

The recommended concentration range for use is 5-20 µg/ml.